

Shiga toxin-producing *Escherichia coli* (STEC) other than the top-7: multiplex PCR assays development and validation for detection and identification of the predominant serogroups in cattle feces.

by

Justin Ludwig

B.S., Kansas State University, 2015

A THESIS

Submitted in partial fulfillment of the requirements for the degree

MASTER OF SCIENCE

Department of Diagnostic Medicine & Pathobiology
College of Veterinary Medicine

KANSAS STATE UNIVERSITY
Manhattan, Kansas

2019

Approved by:

Major Professor
T.G. Nagaraja

Copyright

JUSTIN LUDWIG

2019

Abstract

Shiga toxin-producing *E. coli* (STEC) are major food borne pathogens that cause mild to severe hemorrhagic colitis in humans. Seven serogroups, which include O26, O45, O103, O111, O121, O145 and O157, considered the ‘top-7 STEC’, cause a majority of food borne STEC infections. Cattle are a major reservoir for the STEC. These organisms reside in the hindgut and are shed in the feces. Additionally, cattle harbor a number of other serogroups (n=113) of STEC. A few of the serogroups have been implicated to cause sporadic human infections. Traditionally, serogrouping of *E. coli* is done by agglutination reaction with serogroup-specific antisera. Multiplex PCR (mPCR) assays targeting STEC serogroups, other than the top-7, were designed, optimized, and validated. The assays targeted one or more serogroup-specific genes, *wzx*, *wzy*, *gnd*, *wbdA* of the O-antigen gene cluster. A total of 11 sets of mPCR assays targeting 113 serogroups of STEC were designed and validated with known strains of STEC. The assays were utilized to identify serogroups of *E. coli* strains (n=358), positive for *stx1*, *stx2* (Shiga toxins 1 and 2), and or *eae* (intimin) genes, and negative for the top 7 serogroups, isolated from fecal samples collected from commercial; feedlots. Of the 358 strains from multiple feedlots located in two states, 298 (83.2%) strains belonged to 21 serogroups and 60 strains (16.8%) were unidentifiable. The twenty one serogroups included O168 (29.9%), O109 (17.9%), O131 (8.1%), O2 (7.0%), O104 (4.7%), O171 (4.2%), O74 (3.6%), O8 (2.8%), O136 (0.84%), O178 (0.56%), O102 (0.56%), O175 (0.56%), O98 (0.28%), O118 (0.28%), O113 (0.28%), O96 (0.28%), O76 (0.28%), O35 (0.28%), O20 (0.28%), O160 (0.28%), and O169 (0.28%). The detection of the STEC strains indicated that serogroups, O2, O74, O109, O131, O168, and O171, were the most predominant. A six-plex PCR was designed and validated to detect the prevalence of the top six major non top-7 STEC in cattle feces from feedlot cattle. Among the six

serogroups of the major non-top 7 STEC, serogroups of O2 (48.7%), O109 (88.1%), O168 (67.9%), and O171 (80.4%) were the dominant serogroups in cattle feces collected from commercial feedlots. Fecal samples that contained the six serogroups were also positive for the three major virulence genes, *stx1*, *stx2*, and *eae*. In conclusion, cattle harbor and shed in the feces a number of serogroups of STEC other than the top 7, with O2, O109, O168, and O171 being the predominant STEC. The potential risk of these STEC to cause human infections is not known.

Table of Contents

List of Figures.....	vi
List of Tables.....	vii
Acknowledgements.....	viii
Dedication.....	ix
Chapter 1. An overview of Shiga toxin-producing <i>Escherichia coli</i> (STEC) other than the top-7, in cattle.....	1
Introduction.....	1
Shiga toxin-producing <i>Escherichia coli</i> other than the top-7 in cattle.....	3
Detection methods of STEC other than the top-7.....	5
Conclusion.....	6
Chapter 2. Development and validation of multiplex PCR assays to identify ‘non-top-7’ serogroups of Shiga toxin-producing <i>Escherichia coli</i> known to be present in cattle feces.....	11
Introduction.....	11
Materials and methods.....	12
Primer Design.....	13
PCR assay conditions and optimization.....	13
Validation of PCR assays.....	13
Results.....	14
Discussion.....	14
Chapter 3. Detection of serogroups, other than the top-7, of Shiga toxin gene carrying <i>Escherichia coli</i> (STEC) and prevalence of serogroups O2, O74, O109, O131, O168, and O171 in cattle feces.....	53
Introduction.....	53
Materials and methods.....	54
Non-top 7 STEC strains.....	55
Detection of serogroups.....	55
Development and validation of mPCR targeting O2, O74, O109, O131, O168, and O171.....	56
PCR running conditions.....	56
Validation of the mPCR assay with pure cultures and feces spiked with pure cultures.....	57
Subtyping of the stx genes.....	58
Results.....	59
Discussion.....	60
References.....	68

List of Figures

Figure 1. Multiplex PCR assay (Set 1) to detect serogroups O4, O86, O84, O156, O80, O168, O91, and O109 of Shiga toxin-producing <i>Escherichia coli</i>	15
Figure 2. Multiplex PCR assay (Set 2) to detect serogroups O128, O138, O123, O120, O119, O13, O175, O171, O22, and O5 of Shiga toxin-producing <i>Escherichia coli</i>	16
Figure 3. Multiplex PCR assay (Set 3) to detect serogroups O6, O98, O181, O75, O33, O116, O150, O79, and O25 of Shiga toxin-producing <i>Escherichia coli</i>	17
Figure 4. Multiplex PCR assay (Set 4) to detect serogroups O78, O2, O146, O76, O178, O126, O113, O118, O15, and O147 of Shiga toxin-producing <i>Escherichia coli</i>	18
Figure 5. Multiplex PCR assay (Set 5) to detect serogroups O62/O68, O7, O163, O136, O8, O92, O87, O55, and O20 of Shiga toxin-producing <i>Escherichia coli</i>	19
Figure 6. Multiplex PCR assay (Set 6) to detect serogroups O141, O153, O132, O130, O108, O96, O88, O107/O117, O74, O38, O39, and O115 of Shiga toxin-producing <i>Escherichia coli</i>	20
Figure 7. Multiplex PCR assay (Set 7) to detect serogroups O70, O53, O69, O51, O17/O44/O73/O77/O106, O43, O40, O37, O35, O28/O42 (O20ac), O18 (ab/ac), and O1 of Shiga toxin-producing <i>Escherichia coli</i>	21
Figure 8. Multiplex PCR assay (Set 8) to detect serogroups O139, O125(ab/ac), O124/O164, O90/O127, O102, O105, O85, O82, O81, O148and O140 of Shiga toxin-producing <i>Escherichia coli</i>	22
Figure 9. Multiplex PCR assay (Set 9) to detect serogroups O169, O161, O154, O114, O110, O93, O149, O49, and O21 of Shiga toxin-producing <i>Escherichia coli</i>	23
Figure 10. Multiplex PCR assay (Set 10) to detect serogroups O165, O160, O182, O179, O46/O134, O177, O176, O174, O172, O170, O159, and O152 of Shiga toxin producing <i>Escherichia coli</i>	24
Figure 11. Multiplex PCR assay (Set 11) to detect serogroups O131, O19, O16, O63, O23, O29, O101/O162, O112(ab), O11, O10, and O3 of Shiga toxin producing <i>Escherichia coli</i>	25

List of Tables

Table 1. Shiga toxin-producing <i>E. coli</i> , other than the top-7 reported in cattle feces.....	7
Table 2. Shiga toxin-producing <i>Escherichia coli</i> serogroups identified in cattle.....	26
Table 3. Multiplex PCR (mPCR) assay sets, serogroups and serogroup-specific genes targeted, and size of amplicons.....	27
Table 4. Multiplex PCR assays running conditions for the detection Shiga toxin-producing <i>Escherichia coli</i> (STEC) serogroups, other than the top-7 STEC.....	35
Table 5. Validation of multiplex PCR assays to detect ‘non-top-7’ Shiga toxin producing <i>Escherichia coli</i>	39
Table 6. Serogroup and virulence gene identification of enterohemorrhagic <i>Escherichia coli</i> and Shiga toxin-producing <i>Escherichia coli</i> strains isolated from commercial feedlot cattle feces (2013 study) by end-point PCR....	59
Table 7. Serogroup and virulence gene identification of enterohemorrhagic <i>Escherichia coli</i> and Shiga toxin-producing <i>Escherichia coli</i> strains isolated from commercial feedlot cattle feces (2014 study) by end-point PCR.....	60
Table 8. Target genes, primers sequences, and amplicon sizes to detect six major non top-7 Shiga toxin-producing <i>Escherichia coli</i>	62
Table 9. Fecal spike validation of six major non top-7 Shiga toxin-producing <i>Escherichia coli</i> from cattle feces by end-point PCR method of detection.....	63
Table 10. Detection of six major non top-7 Shiga toxin-producing <i>Escherichia coli</i> and virulence genes, <i>stx1</i> , <i>stx2</i> , and <i>eae</i> , from cattle feces (2013 study) in a commercial feedlot setting by end-point PCR.....	64
Table 11. Prevalence of six major non top-7 Shiga toxin-producing <i>Escherichia coli</i> and virulence genes, <i>stx1</i> , <i>stx2</i> , and <i>eae</i> , from cattle feces (2014 study) in a commercial feedlot setting by end-point PCR.....	65
Table 12. Number and percentage of samples testing positive to multiple major non top-7 Shiga toxin-producing <i>Escherichia coli</i> serogroups (O2, O74, O109, O131, O168, O171) by end-point PCR.....	67

Acknowledgements

The research was supported by the U.S. Department of Agriculture National Institute for Food and Agriculture Grant No. 2012-68003-30155. The author would like to thank Dr. T.G. Nagaraja, Dr. Jianfa Bai, and Dr. Randy Phebus for their many contributions to the thesis research, as well as for their continual demonstration of exemplary mentorship both inside and outside the laboratory. The author would also like to thank Xiaorong Shi and Neil Wallace for their tremendous technical expertise and the numerous graduate and undergraduate students that assisted with the thesis research.

Additionally, the author would like to thank Jim and Patricia Ludwig for their continued encouragement, support, and advice. Many sleepless and stressful nights were beaten by their fortification of belief. Further, the author would not have been able to complete this work without the wonderful aid of their friends and loved ones. Thank you so much to Lance Noll and Kaylen Capps for their contributions and efforts during this work. The author is forever grateful for their collaborations.

Dedication

This work is dedicated to the Ludwig and Behrman Family. Their belief in higher education and continuous learning is just one example of how they allow me to stand on the shoulders of giants.

Chapter 1 - An overview of Shiga toxin-producing *Escherichia coli* (STEC) other than the top-7, in cattle.

Introduction

Escherichia coli, belonging to the *Enterobacteriaceae* family, is a Gram negative, facultatively anaerobic, rod shaped organism. The species is present in various habitats, including gastrointestinal tract of animals and humans, soil and water. Most strains of *E. coli* are non-pathogenic, however, some have become pathogenic because of acquisition of virulence genes via mobile genetic elements. The pathogenic *E. coli* or pathotypes are broadly categorized into those that cause intestinal infections and extra-intestinal infections. The categorization of *E. coli* pathotypes is based on the virulence mechanisms, such as patterns and effects of bacterial attachment onto host cells, effects of attachment on host cells, production of toxins, and invasiveness. The pathotypes that cause intestinal infections, also called diarrheagenic, include at least five major categories: Shiga toxin-producing *Escherichia coli* (STEC) of which enterohemorrhagic *E. coli* (EHEC) is a subset, enteropathogenic *Ea coli* (EPEC), enterotoxigenic *E. coli* (ETEC), and enteroaggregative *E. coli* (EAaggEC) of the pathotypes, Shiga toxin-producing *E. coli* is the most common in terms of prevalence, particularly in the gastrointestinal tract habitat of cattle.

Shiga toxin-producing *E. coli* (STEC) are major zoonotic, foodborne pathogens that are associated with various disease outcomes in humans. Infections caused by STEC can range from mild enteritis to hemorrhagic colitis (HC) to complications of hemolytic uremic syndrome (HUS), particularly in children. The food borne infections are caused by consumption of contaminated beef, produce, or water. Cattle are generally considered major reservoirs for

STEC, however, other ruminant and domestic animal reservoirs have been linked to STEC infections (Gannon et al. 1993). The STEC pathotype possesses a variety of virulence factors that are critical to their adherence to host cells and cytotoxicity. Of the many, Shiga toxins are major virulence factors. There are two types of Shiga toxins, Shiga toxin1 and Shiga toxin-2, encoded by *stx1* and *stx2* genes, respectively. Additionally, there are a number of other virulent factors associated with STEC. The intimin, which is encoded by the *eae* gene, is responsible for adhesion to host cells and in the formation of attaching and effacing lesions in the intestinal epithelium. The possession of this gene distinguishes a subset of STEC called enterohemorrhagic *E. coli* (EHEC; Fratamico et al. 2011, Fagan et al. 1999).

Based on the antigenicity of cell surface structures, *E. coli* is classified into serogroups and serotypes. Serogrouping of *E. coli* is based on the O-antigen of the lipopolysaccharide, which is part of the outer membrane, and is composed of repeating units of one or more sugars or their derivatives. Serotyping includes identification based on the H-antigen, which is the flagellar protein. The importance of this classification system allows for distinguishing strains during outbreaks and for surveillance (DebRoy et al. 2016). The “gold standard” for *E. coli* serotyping based on agglutination reactions, carried out on slides, in tubes, or in 96-well microtiter plates (Debroy et al., 2011).

Seven serogroups of STEC, O26, O45, O103, O111, O121, O145, and O157 (Hara-Kudo et al. 2016), often called ‘top-7’, have been associated with a majority of human infections. The top 7 serogroups have been shown to colonize the hindgut of cattle and are shed in the feces, which is the source of contamination of food products (Tozzi et al., 2003). Among the top 7 serogroups,

the most common serotype is O157:H7. In 1994, the USDA Food Safety and Inspection Service (FSIS) declared *E. coli* O157:H7 as an adulterant in raw ground beef in response to a large outbreak associated with undercooked ground beef. The FSIS began verification testing of ground beef for this pathogen (Wasilenko et al. 2014). In addition to O157, the USDA-FSIS declared the remaining top six serogroups as adulterants in 2012. According to the Centers for Disease Control and Prevention, STEC serogroups O26, O45, O103, O111, O121, and O145 cause a majority of cases of illness due to non-O157 STEC in the US and some of them are major STEC serogroups in other countries (Fratamico et al. 2012). The morbidity and mortality associated with large outbreaks of STEC disease have highlighted the threat these organisms pose to public health.

There have been a total of 184 serogroups of *E. coli* identified, illustrating the vast diversity of *E. coli* (Bettelheim, 2007, Iguchi et al. 2015; DebRoy et al., 2016). A number of serogroups of STEC, other than the top-7, have been identified in cattle (n=108). These serogroups are also harbored in the hindgut and are shed in the feces. According to Ennis et al. (2012) and Monaghan et al. (2011) there have been reports of STEC other than the top-7 linked to human infections. Additionally, Aidar-Ugrinovich et al. (2007) and Gonzalez et al. (2016) have linked STEC serogroups not within the top-7 to causal contaminants of food products and human infections.

Shiga toxin-producing *Escherichia coli* other than the top 7 in cattle.

A total of 113 serogroups of STEC have been identified in feces (Bettelheim, 2007). O1, O2/50, O3, O4, O5, O6, O7, O8, O10, O11, O15, O16, O17/O44/O73/O77/O106, O18 (ab/ac), O19,

O20, O21, O22, O23, O25, O29, O35, O37, O38, O39, O40, O42, O43, O46/O134, O49, O51, O53, O55, O62/O68, O63, O69, O70, O74, O75, O76, O79, O80, O81, O82, O84, O85, O86, O87, O88, O90, O91, O92, O93, O96, O98, O101/O162, O102, O105, O107/O117, O108, O109, O110, O112, O113, O114, O115, O116, O118/O151, O119, O120, O124, O125 (ab/ac), O126, O128, O130, O132, O136, O139, O140, O141, O146, O148, O149, O150, O152, O153, O154, O156, O159, O160, O161, O163, O165, O168, O169, O170, O171, O172, O174, O175, O176, O177, O178, O179, O181, O182 (Table 1). These serogroups are believed to be colonized in the hindgut and are shed in the feces, therefore of concern for food safety(Iguchi et al., 2015, Bettelheim, 2007, Aidar-Ugrinovich et al., 2007, Gonzalez et al., 2016, Guth et al., 2003, Meichtri et al., 2004, Monaghan et al., 2011, Padola et al., 2004)..

Of the 113 STEC serogroups that have been identified in cattle, 85 serogroups have been associated with human illnesses. The serogroups include O1, O2/50, O3, O4, O5, O6, O7, O8, O11, O15, O16, O17/O44/O73/O77/O106, O18 (ab/ac), O20, O21, O22, O23, O25, O38, O39, O40, O43, O46/O134, O53, O55, O62/O68, O69, O70, O74, O75, O76, O79, O80, O81, O82, O84, O85, O86, O87, O88, O91, O92, O93, O96, O98, O101/O162, O102, O105, O107/O117, O109, O110, O112, O113, O114, , O115, O116, O118/O151, O119, O120, O125 (ab/ac), O126, O128, O132, O136, O146, O148, O150, O152, O153, O154, O156, O160, O161, O163, O165, O168, O169, O171, O172, O174, O175, O177, O178, O179, O181 (Iguchi et al., 2015, Bettelheim, 2007, Aidar-Ugrinovich et al., 2007, Gonzalez et al., 2016, Guth et al., 2003, Meichtri et al., 2004, Monaghan et al., 2011, Padola et al., 2004, Hornitzky et al., 2005, Fukushima et al., 2004, Beutin et al., 2007, Zweifel et al., 2005, Irino et al., 2005, Werber et al., 2008). Among the list are specific serogroups that are linked to a higher likelihood to cause

human infections. For example, according to Gonzalez et al., 2016, Meichtri et al., 2004, O8 and O165 have been connected to HUS in children. In addition to O8 and O165, O91, O113, O117, O118, O128 and O146 have a higher frequency of being linked to food contamination and health complications, such as gastroenteritis and hemorrhagic colitis (Blanco et al., 2004). Serogroup O178 has also emerged as a non top-7 STEC that is associated with human illnesses. The O178 serogroup has been frequently isolated in Germany (Beutin et al., 2007; Werber et al., 2008; Slanec et al., 2009) and South America (Cergole-Novella et al., 2007; Oliveira et al., 2008) and connected to sporadic cases of HUS in both regions.

Detection methods of STEC other than the top-7

Agglutination and antisera detection

To date, the traditional method of detection for STEC, other than the top-7, is based on serology. This method determines the O and H antigens by agglutination test carried out in microtiter plates and tubes. The O antisera is produced in a number of laboratories across the world and made available for use. Additionally, strains within a serotype can be isolated based on virulence factors that correspond to the serogroup identified (Hornitzky et al., 2005, Fukushima et al., 2004, Beutin et al., 2007, Zweifel et al., 2005, Irino et al., 2005, Werber et al., 2008, Iguchi et al., 2015, Bettelheim et al., 2007, Aidar-Ugrinovich et al., 2007, Gonzalez et al., 2016, Guth et al., 2003, Meichtri et al., 2004, Monaghan et al., 2011, Padola et al., 2004).

Culture based detection

Culture based detection of STEC, other than top-7, is generally with the use of non-selective media, such as MacConkey agar, and confirmation of the serogroup by agglutination and Shiga

toxin gene by PCR. Hornitzky et al. (2001), using the dual technique (culture and PCR) detected STEC in fecal samples and estimated a prevalence of 18.7%. Fecal samples are cultured onto blood agar supplemented with vancomycin-cefixime-cefsulodin, allowing for the detection of hemolytic colonies, then positive colonies were subjected to DNA extraction and PCR assays

PCR detection

Iguchi et al. (2015) have used 162 primary pairs for the identification of O groups. Of the 162 pairs, 147 pairs were used to identify 147 serogroups based on unique O genes and the other 15 pairs were used to identify 15 groups or strains. Each of these groups contained strains with identical or very similar o antigen genes and included a total of 35 O groups. The authors then used the 162 primer pairs to develop 20 multiplex PCR sets, with each set consisting of 6 to 9 primer pairs of different amplicon sizes. Primer pairs targeted unique sequences specific to O group biosynthesis genes, for example *wzx*, *wzy*, *wzm*, and *wzt*. Iguchi et al. (2015) concluded that PCR-based method is a promising tool for the serogrouping of *E. coli* strains.

Conclusion

Shiga toxin-producing *E. coli* is a major zoonotic foodborne pathogen. The diseases caused by STEC range from mild enteritis to hemorrhagic colitis (HC) and serious complications of hemolytic uremic syndrome (HUS). The complication of HUS is particularly seen in children and is fatal if it is not diagnosed early (Tarr et al., 2005). The top 7 STEC (O26, O45, O103, O111, O121, O145, and O157) are most commonly associated with human infections. Cattle harbor a number of other serogroups of STEC other than the top-7. Although there are several detection strategies for STEC that are the top-7, the detection of other STEC is by random

isolation and confirmation of serogroup by serology. A PCR-based method of detection has been described to identify 147 serogroups, which minclude STEC and non-STEC

Table 1. Shiga toxin-producing *E. coli*, other than the top-7 reported in cattle feces

Animal	Type of sample	Diet	Methodology	Serogroup (Prevalence)	Reference
Cattle	Fecal samples	Not reported	Serology (Agglutination with O-group specific antisera)	O5:H-, O91:H21, O113:H21	Hornitzky et al., 2001.
Beef Steers 14-16 months old	Stool samples and rectal swabs- taken at slaughter level	Not reported	Serology (Agglutination with O-group specific antisera)	O8(16 strains), O113(14 strains), O91(4), O171(3), O174(3), O25(2), O112(2), O2, O11, O128, O143, O146	Meichtri et al., 2004.
Cattle	Fecal samples	Not reported	Serology (Agglutination with O-group specific antisera)	O119, O146, O159, O177, O8	Hornitzky et al., 2005.
Calves, heifers and diarrheagenic animals	Fecal samples of calves, heifers and diarrheagenic animals	Not reported	Serology (Agglutination with O-group specific antisera)	O113, O178, O79 O77, O88, O98	Irino et al., 2005.

cattle, ovines	Isolates recovered from humans, cattle, ovines and food (ground beef)	Not reported	Serology (Agglutination with O-group specific antisera)	O4, O8, O20, O26, O111, O118, O118, O171, O177	Mora et al., 2005.
Bovine,	Strains from human food, environment	Not reported	Serology (Agglutination with O-group specific antisera)	O8, O26, O103, O111, O113 and O128.	Bettelheim et al., 2007.
Steers (80)	Rectal fecal samples	Not reported	Serology (Agglutination with O-group specific antisera)	O2, O49, O110, O132 O156, O172	Diarra et al., 2009.
Cattle	Fecal samples	Not reported	Serology (Agglutination with O-group specific antisera)	O2, O26, O63, O148, O149, O174	Scott et al., 2009.
Cattle	Fecal samples	Not reported	Serology (Agglutination with O-group specific antisera)	O2, O113, O131	Monaghan et al., 2011.

Cattle	Fecal samples (from 12 beef farms, beef suckler and dairy herds)	Not reported	Serology (Agglutination with O-group specific antisera)	O2,O3,O26,O33, O69,O76,O88,O113,O118,O136,O150, O153,O171	Ennis et al., 2012.
Cattle	Fecal samples collected from an abattoir and farms	Not reported	Serology (Agglutination with O-group specific antisera)	O1,O2,O8,O55,O84,O91,O113 (8/39)	Mekata et al., 2014.

Chapter 2 - Development and validation of multiplex PCR assays to identify ‘non-top-7’ serogroups of Shiga toxin-producing *Escherichia coli* known to be present in cattle feces.

Introduction

Escherichia coli, a Gram-negative facultative anaerobic and rod-shaped organism, is a common inhabitant of the gut of animals and humans. The species is categorized into serogroups and serotypes based on O-antigen (somatic) and H-antigen (flagellar). The O-antigen, which is a polysaccharide, is part of the lipopolysaccharide portion of the outer membrane, and the H antigen is the flagellar protein (flagellin) of the bacteria. A total of (n=184) serogroups of *E. coli* have been recognized based on O antigens (Iguchi et al., 2014; Debroy et al., 2015). *E. coli* can be non-pathogenic or pathogenic. The pathogenic strains are categorized into several pathotypes based on the virulence mechanisms. Among the pathotypes, Shiga toxin producing *E. coli* (STEC) appears to be the major one based on the number of serogroups that have been shown to produce Shiga toxins or carry genes (*stx*) that code for Shiga toxins. Shiga toxin-producing *E. coli* are important foodborne pathogens that can cause a variety of diseases, ranging from mild to hemorrhagic enteritis that could lead to complications of hemolytic-uremic syndrome (HUS), particularly in children (Karmali et al., 2004). Specific serogroups and serotypes of STEC are associated with enteritis, hemorrhagic colitis, and hemolytic uremic syndrome. (Stenutz et al., 2006, Knirel et al., 2016). *Escherichia coli* O157 and six other serogroups, O26, O103, O111, O121, and O145, often referred to as the ‘top-7’, are considered as major STEC and have been declared as adulterants in ground beef and non-intact beef products (Fratamico et al., 2014). Cattle are a major reservoir of the top-7 STEC, in which the organisms colonize or reside in the hindgut and are shed in the feces (Hussein et al., 2005). Additionally, as many as 113

serogroups of STEC have been identified to be present feces of cattle (Battelheim, 2007; Aidar-Ugrinovich et al., 2007, Blanco et al., 2004, Gonzalez et al., 2016, Guth et al., 2003, Meichtri et al., 2004, Monaghan et al., 2011, Padola et al., 2004). Some of the non-top-7 STEC have been shown to be associated with sporadic human illnesses (Bettelheim et al., 2007).

Traditionally, serogrouping of *E. coli* is done by agglutination reaction with serogroup-specific antisera, often carried out in a few reference laboratories. A broad serogrouping method based on restriction fragment length polymorphism pattern of amplified)-antigen biosynthesis gene cluster (O-AGC) has been described (Coimbra et al., 2000). In recent years, whole genome sequencing and DNA microarray have been used to identify serogroups of *E. coli* (Norman et al., 2015. Lacher et al., 2014; Liu and Fratamico et al., 2006, Joensen et al., 2015.). A PCR-based assay targeting serogroup-specific genes would be a simpler, low-cost alternative method, readily adaptable to any laboratory, to identify serogroups of *E. coli*. A number of PCR-based assays, end point or real time, have been developed and validated for one or multiple clinically relevant serogroups of *E. coli* (Gonzalez et al., 2016, Monaghan et al., 2011, Iguchi et al., 2014, Blanco et al., 2004; Bai et al., 2012;). Previously, we developed three sets mPCR assays to detect top-7 serogroups and three major virulence genes (Noll et al., 2015, Shridhar et al., 2016) Iguchi et al. (2015) have designed 162 primer pairs to identify 147 individual O-serogroups of *E. coli*, which were then used to create 20 multiplex PCR assays. Because the primary focus of our laboratory is STEC in cattle feces, we wanted to develop PCR based assays to detect 113 serogroups of STEC known to be present in cattle feces (Table 2). Therefore, our aim was to develop and validate several sets of multiplex PCR assays targeting serogroup-specific genes to detect 113 serogroups of STEC.

MATERIALS AND METHODS

Primer design

The nucleotide sequence of the targeted gene for each of the 113 *E. coli* serogroups were obtained from NCBI GenBank, and a search was performed to look for additional corresponding sequences using BLASTN on NCBI's Nucleotide collection (nr/nt) database. The *wzx* gene, which encodes for the O-antigen flippase (Liu et al., 2015) required for O-polysaccharide export, *wzy* gene (Liu et al., 2015), which encodes for the O-antigen polymerase required O antigen biosynthesis, and *wbdA* gene (Kido et al., 2000), which encodes for mannosyltransferase, were used to design primers for serogroup-specific identification. The sequences for each serogroup were aligned using ClustalX version 2.0 (Larkin et al., 2007). Primers were designed to amplify the targets with distinct amplicon sizes that can be differentiated on a capillary gel electrophoresis equipment. The primer sequences for the 113 serogroups are shown in Table 3.

PCR assay conditions and optimization

The PCR protocol for specific gene target included an initial denaturation at 94 °C for 5 min, followed by 25 or 30 cycles of denaturation at 94° C for 30 s, annealing for 30 s at 59-68°C, extension for 75 s at 68°C and a final step of extension at 68°C for 7 min. The number of PCR cycles and annealing temperatures varied based on optimization for each set (Table 4). All other conditions were same for all 11 sets of assays (Table 5).

Validation of PCR assays

Strains belonging to 113 serogroups of STEC, other than the top-7 (O26, O45, O103, O111, O121, O145, and O157), were used in the study (Table 3). The strains were obtained from our

culture collection (n=104), *E. coli* Reference Center at Pennsylvania State University (n=231), Michigan State University (n=5), University of Nebraska (n=5), and Food and Drug Administration (n=36). Strains, stored in protect beads at -80° C, were streaked onto blood agar plates (Remel, Lenexa, KS) and incubated at 37 C. Colonies were picked and suspended in 1 ml of distilled water, boiled for 10 min, centrifuged at 9,300 x g for 5 minutes and the supernatant was used for PCR assays. Each set of assays was validated by testing with known strains of STEC (n= 381) belonging to 113 serogroups (Table 5).

Results

A total of 11 mPCR assays, each targeting 8 to 11 serogroups, targeting a total of 113 serogroups of STEC were designed. Each set contained primer pair sequences that generated amplicons of different sizes that can be readily differentiated in a capillary electrophoretic system (Table 3); Figures 1 to 11). The running conditions ranged from 25 to 30 cycles and annealing temperatures varied from 59°C to 66°C (Table 4). Initial optimization and validation of the assay were done with one strain and subsequently validated with additional strains (Table 4).

Discussion

Traditionally, serogrouping and serotyping *E. coli* involved agglutination reactions with specific antisera raised against O and H antigens (DebRoy et al., 2011). This method is labor-intensive and in some instances, specificity is questionable because of cross-reactions with other serogroups (Ennis et al., 2012). The development of PCR assays would be a simpler and effective tool to identify serogroups. More importantly, the PCR assay is readily adaptable to any laboratory. Many PCR assays have been developed and validated, generally targeting 1 to 7

serogroups, often with major virulence genes (Shiga toxins 1 and 2, intimin, and enterohemolysin (Noll et al. 2015). The colonization and fecal prevalence of STEC is widespread among cattle, indicating the diverse possibility of contamination (Blanco et al., 2004). Additionally, the isolation and of non-top-7 STEC can be challenging because the lack of immunomagnetic beads and selective media differentiation (Hornitzky et al., 2001). STEC serogroups other than the top-7 have been previously identified in other studies. Serogroups that have been identified from cattle comprise of 108 serogroups. Of the 108 STEC serogroups that have been identified in cattle, 85 serogroups have been associated with human illness. (Iguchi et al., 2015, Bettelheim et al., 2007, Aidar-Ugrinovich et al., 2007, Gonzalez et al., 2016, Guth et al., 2003, Meichtri et al., 2004, Monaghan et al., 2011, Padola et al., 2004). Serogroups O8 and O165 have been found in association with HUS in other countries, including Brazil and Argentina (Gonzalez et al., 2016, Meichtri et al., 2004). Additionally, humans have a higher likelihood of becoming exposed to O91, O113, O117, O118, O128 and O146 because of their higher prevalence in animals and subsequent contamination of food and water. These serogroups are frequently associated with severe illnesses (Blanco et al., 2004). The design of multiplex PCR assays allows for a higher throughput and increased specificity compared to other methods such as agglutination reactions with antisera, or culture method approach. (Irino et al., 2005, Zweifel et al., 2005).

Figure 1. Multiplex PCR assay (Set 1) to detect serogroups O4, O86, O84, O156, O80, O168, O91, and O109 of Shiga toxin-producing *Escherichia coli*.




Figure 2. Multiplex PCR assay (Set 2) to detect serogroups O128, O138, O123, O120, O119, O13, O175, O171, O22, and O5 of Shiga toxin-producing *Escherichia coli*.




Figure 3. Multiplex PCR assay (Set 3) to detect serogroups O6, O98, O181, O75, O33, O116, O150, O79, and O25 of Shiga toxin-producing *Escherichia coli*.




Figure 4. Multiplex PCR assay (Set 4) to detect serogroups O78, O2, O146, O76, O178, O126, O113, O118, O15, and O147 of Shiga toxin-producing *Escherichia coli*.




Figure 5. Multiplex PCR assay (Set 5) to detect serogroups O62/O68, O7, O163, O136, O8, O92, O87, O55, and O20 of Shiga toxin-producing *Escherichia coli*.




Figure 6. Multiplex PCR assay (Set 6) to detect serogroups O141, O153, O132, O130, O108, O96, O88, O107/O117, O74, O38, O39, and O115 of Shiga toxin-producing *Escherichia coli*.




Figure 7. Multiplex PCR assay (Set 7) to detect serogroups O70, O53, O69, O51, O17/O44/O73/O77/O106, O43, O40, O37, O35, O28/O42 (O20ac), O18 (ab/ac), and O1 of Shiga toxin-producing *Escherichia coli*.




Figure 8. Multiplex PCR assay (Set 8) to detect serogroups O139, O125(ab/ac), O124/O164, O90/O127, O102, O105, O85, O82, O81, O148and O140 of Shiga toxin-producing *Escherichia coli*.




Figure 9. Multiplex PCR assay (Set 9) to detect serogroups O169, O161, O154, O114, O110, O93, O149, O49, and O21 of Shiga toxin-producing *Escherichia coli*.




Figure 10. Multiplex PCR assay (Set 10) to detect serogroups O165, O160, O182, O179, O46/O134, O177, O176, O174, O172, O170, O159, and O152 of Shiga toxin-producing *Escherichia coli*.




Figure 11. Multiplex PCR assay (Set 11) to detect serogroups O131, O19, O16, O63, O23, O29, O101/O162, O112(ab), O11, O10, and O3 of Shiga toxin-producing *Escherichia coli*.




Table 2. Shiga toxin-producing *Escherichia coli* serogroups identified in cattle feces.

O1	O2/50	O3	O4	O5	O6	O7	O8	O10	O11
		O17/O44/							
O15	O16	O73/O77/	O18	O19	O20	O21	O22	O23	O25
		O106							
O29	O35	O37	O38	O39	O40	O42	O43	O46	O49
O51	O53	O55	O62	O63	O69	O70	O74	O75	O76
O79	O80	O81	O82	O84	O85	O86	O87	O88	O90
O91	O92	O93	O96	O98	O101	O102	O105	O107	O108
O109	O110	O112	O113	O114	O115	O116	O118	O119	O120
O124	O125	O126	O128	O130	O131	O132	O136	O139	O140
O141	O146	O147	O148	O149	O150	O152	O153	O154	O156
O159	O160	O161	O163	O165	O168	O169	O170	O171	O172
O174	O175	O176	O177	O178	O179	O181	O182		

Table 3. Multiplex PCR (mPCR) assay sets, serogroups and serogroup-specific genes targeted, and size of amplicons

mPCR set	Serogroup	Genes targetd	Primers sequence	Amplicon size (bp)
1	O4	wzx	F: TGAAACAGCAGTGCCTGCATTCTC R: GATTGCCGTGCCAATTATGCCGAA	832
	O80	wzx	F: GAAGCATGGCTTCTAGGGGG R: TGACAAAGGTAGCCACGGAA	406
	O84	wzx	F: TCAGCGTTCCAAGAAGCACT R: TGGTGTGCACTTATACATCCGA	501
	O86	wzy	F: ATTTGAGGCTGACGCGTATGGACT R: AGCAACACTTCCAATGATCCACCC	562
	O91	wzy	F: CGCATTTAAGGACTGGCTGT R: GTAGCAGATATGCCGACCGT	277
	O109	wzx	F: TCTCTCTCGACATAACCGCGCTT R: ACCGTAGCCCAAAGAGCCACA	204
	O156	wzx	F: TGCTCATGCGTTTAAAATGG R: CAACAAGAACAGGCATCAGG	452
	O168	wzx	F: TGTCGACTTTGGGAAATGTGG R: CTGCAGAGGCCAATTCAGGT	336
	O5	wzy	F: AGCCCAAGAGGCCGCGAGTATTT R: TCCAACCACCATTCTCCGCCAT	176
	O13/O129/ O135	wzx	F: CGGGAGAGCAGTGTTTCCAA R: TATTGCGGCACCCAGTAACC	364
2	O22	wzy	F: TTTACTGGCTGCTGCTAGTGCT R: TCATCTCCACCACGAGTCGAAAGT	246
	O119	wzx	F: CTGGGGCAATCTGCTTTCCT R: CCCAAGGTATTGTTGCCCT	421

	O120	wzx	F: CTGGTTTTGTTGTTGCATTGCT R: CCAGTTGGTGCCAACCAAAG	535
	O123/O186	wzx	F: ACAATTAGGGCCTGGTGCAT R: TGTGCTAGCGCTAAAGGACT	619
	O128	wzx	F: GCCATTACGACGTTGATGACT R: TGCAACCCCAATAGCAAAAGC	768
	O138	wzx	F: GCAGCAATGCCTGCTGTTTT R: AGCGTATGCAACCCCAATGA	696
	O171	wzx	F: TGCTCAAGTGGCATGCAGAT R: TGCAACCTGATATCCAGCAGT	281
	O175	wzx	F: TGTTTAATTGTTCCCTCCGCTCA R: GCAGCCCAACCAAACCTTAGC	343
3	O6	wzy	F: GGATGACGATGTGATTTTGGCTAAC R: TCTGGGTTTGCTGTGTATGAGGC	783
	O25	wzy	F: AGAGATCCGTCTTTTATTTGTTTCGC R: GTTCTGGATACCTAACGCAATACCC	230
	O33	wzx	F: GGACCTCTGTTTGTGGCTGT R: GGTGTGCCTATCGCATAACC	413
	O75	wzy	F: GAGATATACATGGGGAGGTAGGCT R: ACCCGATAATCATATTCTTCCCAAC	511
	O79	wzx	F: GCTATAGTGCACCAGGATTGT R: GATTGCCGTCTGCCCTAAT	266
	O98	wzx	F: GGGCGTATTGAGGTTCTTGT R: CAACTGAAATGGAGCAGCAAAT	675
	O116	wzx	F: CTTTTGCGTTGTGCCTCGAA R: ATACTGCCCTACGTTTGCGG	355
	O150	wzx	F: TGC GTTCACTTGCTGGTTTG R: ATGAGTGCAGGCACTTGAA	313
	O181	wzx	F: TGGAGTAACGAAATACACCGCT R: AGATTGCCAATAACCAGAAGCA	595
4	O2/O50	wzx	F: TGGCCTTGTTTCGATATACTGCGGA	819

			R: TCACGAGCTGAGCGAAACTGTTCA	
O15	wzx	F: CCAAATTGGTTTTTGCAAGG	288	
		R: TGCACTAGTCGCAAGTGTTGA		
O76	wzx	F: CATATGCAGATTGAAGGTAG	533	
		R: GAAAGCCATAAAGTGCC		
O78	wzx	F: GGTATGGGTTTGGTGGTA	992	
		R: AGAATCACAACTCTCGGCA		
O113	wzy	F: GCATGTATGATGCATAGCTTCGCC	419	
		R: TGATATCGTTTCGCTAACCACCCA		
O118/O151	wzx	F: GTGGGAGTCTGAATCAAGTTGCGA	344	
		R: AGCAACCTTACCCAATCCTAAGGG		
O126	wzy	F: CGCATTAATGGACCTGATAAAGCATCG	465	
		R: ACTAGCGCACATATCGTTAGCACG		
O146	wzx	F: AGGGTGACCATCAACACACTTGGA	640	
		R: AGTTCAATACTGTCGCAGCTCCTC		
O147	wzx	F: GCTTGATTGGCGGTAGTGTAT	230	
		R: TGGTAATCCAGCCACAAAAGAA		
O178	wzx	F: CTGTCCGTAGTGAGGTTGGC	495	
		R: ACCTCCAGATCGGTCCTTAATC		
5	O7	wzx	F: GCCAGGCAAGGATTATTATGG	753
		R: TGTCATTGCAGGTACGCTAGA		
	O8	<i>Orf469</i>	F: CCAGAGGCATAATCAGAAATAACAG	448
		R: GCAGAGTTAGTCAACAAAAGGTCAG		
	O20	wzx	F: GAGCAGCGGAATACTTTCCA	204
		R: TCCAAGAGTCCTGATGCAAAT		
	O55	wzx	F: TGCCAATACGTAATACCCAAAA	262
		R: GTGGCCACAGGCAATCTTAT		
	O62/O68	wzx	F: TGAAATGTTTGAGCGTGAATTT	906
		R: CATATCGCCGAATGAGTAACC		
	O87	wzx	F: TGGGTTTTTGTC AAGCATCA	306
		R: TCGCTTTCATTTTCTCCATTC		

	O92	<i>gnd</i>	F: TGC GAAGTGGTCTGAATCTG R: TTCCCAAGTACAATGCACGA	375
	O136	<i>wzx</i>	F: GCCACGTACAAAAATTTAGCC R: CCCCTTCTTGATCCATTTGA	528
	O163	<i>wzy</i>	F: TTGCTGATGATGAGGTTGTGA R: ATGTTCCACTCCCTGCAAGT	596
6	O38	<i>wzx</i>	F: AAGAGGCTGGGGCATTTAGT R: TGCGCCGATAAATCCATATT	253
	O39	<i>wzx</i>	F: CACAAGCATTTCAGCCAAAG R: TGAGAAAATCACCAGCCTGC	201
	O74	<i>wzx</i>	F: CTGGTCAATGGCAAGCTGTA R: ATGCAAAAATCCAAGCCAAT	303
	O88	<i>wzx</i>	F: TTCTCTCCCTTCGTTGGCTA R: TTCCCACACCAGCATTAACA	394
	O96	<i>wzy</i>	F: GGGCTGTTTACGAGAGCATT R: GCACTACTAGGATCATCCGCA	457
	O107/O117	<i>wzx</i>	F: TTTGATGTCGTTGCTCCGTA R: TGCACCAATTGTCGGTGA	357
	O108	<i>wzx</i>	F: GGTGGTCGTGAGCAACTAGG R: TTATACGCGAGATGCTTTGC	515
	O115	<i>wzx</i>	F: GTGAATTTGCTCGCCTCTCT R: TGCTACCAACCATTATGCG	158
	O130	<i>wzx</i>	F: GAGGGCTAATTGCATCCGTA R: GCAACAAATGAACGCATGAC	567
	O132	<i>wzx</i>	F: TGGCAATCCCTCTGATCTTC R: ATCAAACATAACCCGCCTGA	652
	O141	<i>wzx</i>	F: GTGTATTTAACATTTGCTCAAGCC R: ACGCTGAAGGAATCGTCAAC	880
	O153	<i>wzx</i>	F: TTCGTCCTATGCGCAGTATTC R: ATTGCCTGAACACGGAAGAT	741
7	O1	<i>wzx</i>	F: CTTTTAGTGCAAGTCGTGCG	152

			R: AGCTTGCACAGGCAAAAGAT	
O18	wzx	F: CGAGAAGTCGCAATTGAAAA	199	
		R: TGGATGCTGCGAGAATTTTA		
O28	wzx	F: CTCAAGTCATTGGCGCATT	255	
		R: GGCTGACTGGGGTCGTTAT		
O35	wzx	F: TGGCATAGCGGTGTTCTATG	305	
		R: CCCAAGAAAGTGGAACAA		
O37	wzx	F: TAGCCATATGCCTTGTCGTG	353	
		R: TGAATCGCTGGTAACGAACA		
O40	wzx	F: CCTGAAGAGGCTGGGAGAAT	396	
		R: AAGAAAACCTCATTGGCACACC		
O43	wzx	F: TTGAATCGAAGGCTTTTGC	445	
		R: AGGTATCACACCATGCACGA		
O44/17/73/7 7/106	wzx	F: GGAGAAGATGTAGGCGGAGA	500	
		R: AGCACATCAAAAGCTGACGTT		
O51	wzy	F: TTGGTCGCCTTATTTATTGCTT	566	
		R: AACGTAATGAAAAACCAGGAACA		
O53	wzx	F: CGACCTCCTTTATTGGTTGC	735	
		R: TCTTTCCCATAATTCAATCCAA		
O69	wzx	F: ATCATTTGCAGGGCTTATGG	649	
		R: CTCCTTTCATCCCGAGAAT		
O70	wzx	F: ATACCGGTCCATCAGGAATG	863	
		R: AACATCCCGCGAAAAATAAA		
8	O81	wzx	F: AGATCGCAGGCTTCGTTG	248
		R: GGGCAGGACCTTCTTAAACA		
	O82	wzx	F: TCGGGATATCTGGATTTGGA	301
		R: AGCCAATTCATTTGGAATGTTT		
	O85	wzx	F: TGGATACTGCAGGTGAGTGG	353
		R: AATGCTAAGTACAGCGCAACC		
	O90	wzx	F: CAGTTCGGCTGGTTGGTT	498

			R: CCCAGTTGTGCAACAATAACA	
	O102	wzx	F: GCTCGTGAATTGGGAGCA	453
			R: TGCACATGTCAGGGTTGG	
	O105	wzx	F: TCAAAGATGAGAGTGGGAATAGC	407
			R: CGGCTCTTCATCATTGCAT	
	O124/O164	wzx	F: CCGCGATGAATGATTCTGTA	570
			R: TCACTCCCTGAACTCTGCAC	
	O125	wzx	F: GAGAGGAGCATTGATTGCTGT	652
			R: TTGAGCATTTTCGACGTAGC	
	O139	wzx	F: GGGATTAGGCGGCTTTTAAAT	859
			R: TACTTGCCCCAAAAGAGACG	
	O140	wzx	F: AAACGCAATAAGTGGGTG	155
			R: GCAATTGGACCAGAAACAAA	
	O148	wzx	F: TCGCTAAACTGACAGGTGTGAT	201
			R: ACGGCCACGCTTCTTTTA	
9	O21	wzx	F: GCTCAATATGAGTGAGGCTGGT	145
			R: AAAAACGCCTCGTGCTTTTA	
	O49	wzx	F: TTTACCGCGAATAACGGTTC	197
			R: ACAGGGGCAATACGTTTCATC	
	O93	wzx	F: GGCAGGAGAATATGCAAGGA	299
			R: TAAATGCCTCCCAAACACCA	
	O110	wzx	F: TGTTGCAAGCATGAACTTGA	346
			R: TAAGCCATTTTCGCCAAATC	
	O114	wzx	F: TGCTACAAGTGGTGGATTGC	396
			R: CCCATACAACCATCGCAAAT	
	O149	wzx	F: AGCGGTGCAAAGTTAATTCC	253
			R: CACCAAACACATTCTGCGTAA	
	O154	wzx	F: TTGGTTTTTCCAAGGGATACA	499
			R: CAAAGCAAATACCACCATCG	
	O161	wzx	F: GGATATCTGGTTGGCGACAT	646
			R: AGCTTAAAACATGAGGGAGCTG	

	O169	wzx	F: GAAGGAATGGGATTTGAAAAGA R: CGCTTTAACAATTGCTTTCG	865
10	O46/O134	wzx	F: TCAGGTGCGCCATTATTTTT R: TCACCCCCATAATAACCATCTT	455
	O152	wzx	F: TCTGTCGCTATTGTCTTTTCACA R: TGCCGTTAACGTACCAAAAA	150
	O159	wzx	F: AAATCCATTGGTGGAGTAGGAA R: CTTCCGCAAATAAGGTCCAA	202
	O160	wzx	F: TGGCGTTCTAGGTCTTGCTAT R: TGCAGCCCCCTCTGATTCTAT	655
	O165	wzx	F: TGCATTGCTTTATTTGAGC R: TCACGCTTTAACGCATACAGA	735
	O170	wzx	F: AATAGAATCGCCTGGGGATAG R: CAGTGCTTAAAGGATGGTTGG	233
	O172	wzx	F: CGGGCATGTTTGTTTTTGAT R: AGCATCCCACCTCTCACAAT	278
	O174	wzx	F: GGAGATAAAGCTGCAGGTGAG R: CCAGTAAGCGGGCCTAAAAG	317
	O176	wzx	F: TTTTATAGATTTGTTGATCCATTGGT R: CATAAGCGGCACTGCAATTA	356
	O177	wzx	F: GGGGTATACACGCTCGCTAT R: CTTACAGCAGGCAACGCAAC	395
	O179	wzx	F: ATGCCATTGCCGATGATATG R: ACCTCCACCCCAGCATAACT	505
	O182	wzx	F: TTAGTTAATGCAGCCGTCACA R: GCCGATACTCCCATAAAAAGC	566
11	O3	wzx	F: CGCTTAATTCGACCCAAATC R: AGTCGTTCTTGAGCCATCTGTA	145
	O10	wzx	F: AGGAAAACATTTGCCGGAAT R: TATGGAGCGGTGATGAAACA	187
	O11	wzx	F: CCATTGGTCACACTGCCTTA	225

		R: CAAACCAAGAGATGCTGCAA	
O16	wzx	F: TTTTGGGCGGTGATATTTTC	505
		R: CCGCATTAATAATCCCTACG	
O19	wzx	F: GCGGAATATCTCGTGCAGTT	574
		R: GCGAAATGTGAGACAGCAAA	
O23	wzx	F: TTTATCTTGGGCGGCTTAAA	403
		R: CCCACCCAAGCATTTACTTT	
O29	wzx	F: TTATTGGGTTGTCGCCATTT	348
		R: CGGTCAACTTGTGAAAGCAA	
O63	wzx	F: GGGCAGATGCAAGGATTTAC	455
		R: AAGACTCCTCCCCCAATCAG	
O101/O162	wzm	F: TGGACCACATTAGGTTGGTTAGAG	309
		R: AAGCATCATAAAATTTGCGCCATAG	
O112	wzx	F: ACAACGTGAAATTTTTGTGTGG	270
		R: GCACGTATGCATTCTCAAGC	
O131	wzx	F: GTGATTTCTGGGGCAACATT	655
		R: AAGCCTGCCCTAAACAAAGC	

Table 4. Multiplex PCR assays running conditions for the detection Shiga toxin-producing *Escherichia coli* (STEC) serogroups, other than the top-7 STEC.

Assays	Number of O groups	PCR cycles	Annealing temperature	O groups
Set-1	8	25	65°C	O4, O80, O84, O86, O91, O109, O156, O168
Set-2	10	30	65°C	O5, O13/O129/O135, O22, O119, O120, O123/O186, O128, O138, O171, O175
Set-3	9	30	64°C	O6, O25, O33, O75, O79, O98, O116, O150, O181
Set-4	10	30	63°C	O2/O50, O15, O76, O78, O113, O118/O151, O126, O146, O147, O178
Set-5	9	30	61°C	O7, O8, O20, O55, O62/O68, O87, O92, O136, O163
Set-6	12	30	66°C	O38, O39, O74, O88, O96, O107/O117, O108, O115, O130, O132, O141, O153
Set-7	12	30	63°C	O1,O17/O44/O73/O77/ O106, O18 (O18ab/O18ac), O28(O28ac/O42), O35, O37, O40, O43, O51, O53, O69, O70
Set-8	11	25	68°C	O81, O82, O85, O90/O127, O102, O105, O124/O164, O125 (O125ab/O125ac), O139, O140, O148
Set-9	9	25	63°C	O21, O49, O93, O110, O114, O149, O154, O161, O169
Set-10	12	25	59°C	O46/O134, O152, O159, O160, O165, O170, O172, O174, O176, O177, O179, O182

Set-11	11	30	62°C	O3, O10, O11, O16, O19 O23, O29, O63, O101/O162, O112, O131
Total	113			

Table 5. Validation of multiplex PCR assays to detect ‘non-top-7’ Shiga toxin producing *Escherichia coli*

mPCR set	Serogroup	Source of strains					No. of strains positive/No. of strains tested
		Kansas State University	Pennsylvania State University	Michigan State University	University of Nebraska	Food and Drug Administration	
1	O4	6282	O4 std., 14.0484, 14.0486, 14.1232				5/5
	O80		O80 std., 13.0406, 12.3395, 9.1802, 9.1803				5/5
	O84	10340-1, 12824-1					2/2
	O86		O86 std., 14.0607, 14.0615, 14.0747, 14.0748				5/5
	O91	4162, sPRH-561					2/2
	O109	12662-2, 15150-1, 15166-1	13.0166		E10F-6		5/5
	O156		O156 std., 15.0353,				4/4

		15.0354, 15.0355	
	O168	O168 std., 15.0133, 15.0181, 15.0185, 14.0047	5/5
2	O5	14.0590, 14.1660	2/2
	O13/O129/ O135	13.0102, 12.2598	2/2
	O22	1245, 4590, 4570, 1027-4	5/5
	O119	15.1553	1/1
	O120	O120 std., 13.0021, 13.0022, 13.0023, 13.0024	5/5
	O123/O186	O123 std., 8.0269, 7.0571, 5.0769, 15.0862	5/5
	O128	O128 std., 13.0062, 13.0027, 13.1024, 14.0696	5/5
	O138	O138 std., 13.1113	2/2

	O171	1044-1, 12772		2/2
	O175		O175 std., 9.1703, 10.0047, 10.0048, 10.0472	5/5
3	O6	5750-1, 16857, 8276-1		3/3
	O25		O25 std, 3.1421, 14.1138, 14.1119, 14.0200, 14.0052	6/6
	O33		O33 std., 13.0050, 13.1199, 15.1254, 16.0504	5/5
	O75		O75 std., 11.1836, 12.0940, 12.0944, 12.1624	5/5
	O79		O79 std., 5.2457	2/2
	O98	sPRH-571	13.1267, 15.0710	3/3
	O116	1238-1, 1582-1, 99- 3708, 3702-1,		7/7

7712-1,
7716,
7752-2

O150	12.3283, 12.3280, 12.3284, 3.4896, 3.4897	5/5
------	---	-----

O181	O181 std., 11.0228, 11.0392, 12.0697	4/4
------	---	-----

4	O2/O50	99-18866, 99-18858	14.1652	3/3
----------	--------	-----------------------	---------	-----

O15	2794, 18925-1, S918, s915, N5789, N15018	O13 std	TW1681	8/8
-----	--	---------	--------	-----

O76	O76 std., 14.0706, 14.0760, 14.0940, 14.0986	5/5
-----	--	-----

O78	ATCC 35401, ATCC 43896	O78 std., 10.0182	4/4
-----	---------------------------------	----------------------	-----

O113	1598-2, 1862-1, 3517	3/3
------	----------------------------	-----

O118/O151	KDHE 14, KDHE 52	13.1018	3/3
-----------	---------------------	---------	-----

	O126	99-17409	16.0457, 16.0511, 13.0069	4/4
	O146		7.2907, 7.2663, 7.2803, 7.3159, 14.0976	5/5
	O147	13466-1	TW7464	2/2
	O178		O178 std., 14.0800, 14.0828, 14.0829, 14.0908	5/5
5	O7		16.0405, 14.0598, 14.0551	3/3
	O8	2089-2, 4202, 4950-2, sPRH-559, 5702-1, 6102-2, 6504-1, 6846-2, 6849-1, 15484- 1,4544-1, 4964-1, 18370-1, 4542-1, 6412-1, 6416-1,	86.1006	19/19

		5344-1, 5799	
	O20	5481	1/1
	O55	5906, ATCC 12014	TW04062, TW08260 4/4
	O62/O68	O68 standard, 1.2557, 4.0175, 4.2378	4/4
	O87	15.0108, 14.0979, 12.3233	3/3
	O92	13.0124	1/1
	O136	4709-1, 10314-1, 10320-1	14.0889 4/4
	O163	9388-1, 13802-1, 18917-1	3/3
6	O38	748-1, 13950-1, 2072-3	3/3
	O39	13.1182, 15.1383, 16.1097	3/3
	O74	4558-1, 8608, 13472-1, 1229-1	4/4
	O88	1235-1, 2076-2, 8652-1,	5/5

	4560-1, 1240-2			
O96	497, 3712- 1, 18862-1, 3714-1			4/4
O107/O117	3536-3, 9966-1, 3534-1			3/3
O108	9924-1			1/1
O115		13.1261		1/1
O130	492-1, 1038-2, 1239-2, 3270-1			4/4
O132	2067-1, sPRH-568			2/2
O141		13.1660, 14.0312, 14.0314		3/3
O153	1932, 9916-1			2/2
7	O1	14.0223	6227	2/2
	O17/O44/O 73/O77/O10 6	8.2525, 8.2526, 14.0994, 14.0989, 14.0398	3829	6/6
	O18	12.0213, 13.0442, 13.1139	6853	4/4
	O28	4.0349, 5.0190	6404-2	3/3
	O35	1.2208, 6.089	6709	3/3

	O37	5.0628, 7.3209	6594	3/3
	O40	12.0572, 13.0486		2/2
	O43	12.0257, 15.1752, 16.0015	6861	4/4
	O51	11.1861, 12.2022, 15.0773	3724	4/4
	O53	12.0583	6503	2/2
	O69	14.1063, 14.1064	3560	3/3
	O70	95.1266, 96.1774	6388	3/3
8	O81	14.1380, 14.1397	6754	3/3
	O82	14.0846, 15.0429	5541	3/3
	O85	15.1549, 16.0447	3620	3/3
	O90/O127	12.3598	TW4549	2/2
	O102	12.3542, 12.3624, 14.0030	3558	4/4
	O105	9.0695, 9.0698	6810	3/3
	O124/O164	15.0880, 15.0790		2/2
	O125	13.0499, 13.1058	6517	3/3
	O139	13.0785, 14.0081	6672	3/3
	O140	12.066	5737	2/2

	O148			12.3461, 13.1357	2/2
9	O21		9.0069, 12.3518	6508	3/3
	O49		13.0005	5598	2/2
	O93		87.1672	6613	2/2
	O110		12.2044	3744	2/2
	O114		12.0637, 13.1247	6583	3/3
	O149		12.2184, 13.1187	6495	3/3
	O154		13.0130, 14.0265	6742	3/3
	O161			6341	1/1
	O169		12.0689	6992	2/2
10	O46/O134	1223-4	12.1020, 15.0242, 16.0631	3606	5/5
	O152		4.2414, 7.0560		2/2
	O159	1554-1	12.1465		2/2
	O160		15.1575, 14.0525		2/2
	O165	7050-4	13.0107, 14.1095	6757	4/4
	O170		13.0772, 14.1230	6833	3/3
	O172	9042-1, 11154-1			2/2
	O174		14.153	6129	2/2
	O176		14.0193	6996	2/2
	O177		14.1088, 14.1094	E28F-5, 89H-4, 58H-1	5/5

	O179		12.0472, 12.0696, 11.1855, 11.1606		4/4
	O182	sPRH-569	13.151, 14.1296	5745	4/4
11	O3		14.1045, 14.1048	6291	3/3
	O10		7.1714	6888	2/2
	O11		13.0841, 13.0861	6812	3/3
	O16			6734	1/1
	O19		14.0988, 15.0093, 16.054		3/3
	O23		1.2315, 3.2841	6556	3/3
	O29	10720-1	1.1116	3608	3/3
	O63		3.2810, 3. 4555	6441	3/3
	O101/O162			E907F-2	1/1
	O112		14.1406, 14.1409	5127	3/3
	O131		12.3205, 13.0803	6901	3/3

Chapter 3 - Detection of Serogroups, Other than the Top-7, of Shiga Toxin-Gene Carrying *Escherichia coli* (STEC) and Prevalence of Serogroups O2, O74, O109, O131, O168, and O171 in Cattle Feces.

Introduction

Shiga toxin-producing *Escherichia coli* (STEC) are food borne pathogens that induce mild to severe bloody diarrhea (hemorrhagic colitis), which could lead to complications of hemolytic uremic syndrome (HUS), particularly in children. Cattle are major reservoirs of STEC, which reside in the hindgut and are shed in feces. Shiga toxins (Stx) are primary virulence factors associated with STEC induced illness in humans. Individual differences in amino acid sequences of Stx proteins, Stx1 and Stx2, or nucleotide sequences in *stx* genes, *stx1* and *stx2*, as well as differences in cytotoxicity make it possible to identify Stx (*stx*) subtypes (Strockbine et al., 1986; Weinstein et al., 1988; Basu et al., 2015). Subtypes of *stx1* and *stx2* genes include *stx1a*, *stx1c*, and *stx1d*, and *stx2a*, *stx2b*, *stx2c*, *stx2d*, *stx2e*, *stx2f*, and *stx2g*, respectively. Outcomes of human illness are influenced by subtype. Shiga toxin 2 is more commonly associated with HUS than Stx1 (Ethelberg et al., 20014) and of the Shiga toxin subtypes, *stx2a*, and *stx2c* are more commonly associated with HUS (Persson et al., 2007). Intimin, encoded by *eae*, characteristic of enterohemorrhagic *E. coli* (EHEC), is responsible for adhesion of organisms onto intestinal epithelial cells causing attaching and effacing lesions, which is commonly seen in cases of severe human illnesses. The presence of *eae* and *stx2* together result in higher risk for severe infections. However, Stx2 can cause severe illness without the the intimin gene as was seen in the O104:H4 outbreak in Germany in 2011 (Karch et al., 2012). Additionally, *E. coli* hemolysin, also called α -hemolysin (Cavalieri et al., 1984), EhxA, is also an important virulence factor that belongs to

the repeats-in-toxin (RTX) family and lyses red blood cells in the host, which could contribute to severity of illness induced by STEC or EHEC (Welch et al., 2005).

Escherichia coli O157 and six other serogroups, O26, O103, O111, O121, and O145, often referred to as the ‘top-7’, are considered as major STEC and have been declared as adulterants in ground beef and non-intact beef products. Additionally, as many as 108 serogroups of STEC have been identified to be present feces of cattle (Battelheim, 2007; Aidar-Ugrinovich et al., 2007, Blanco et al., 2004; Gonzalez et al., 2016; Guth et al., 2003; Meichtri et al., 2004; Monaghan et al., 2011; Padola et al., 2004). Some of the non-top-7 STEC have been shown to be associated with sporadic human illnesses (Bettelheim et al., 2007). Not much is known about the prevalence of the non-top-7 STEC in cattle feces. In our previous studies (Dewsbury et al., 2015; Cull et al., 2017) that were designed to study the prevalence of the top-7 STEC, a number of *E. coli* strains that were positive for Shiga toxin gene but negative to the top-7 serogroups were obtained. Therefore, our primary objective was to use previously developed multiplex PCR (mPCR) assays to identify the ‘non-top-7 STEC’ strains isolated from feedlot cattle feces and to determine subtypes of *stx* genes.

Because the detection of the STEC strains from the two feedlot studies indicated that serogroups, O2, O74, O109, O131, O168, and O171, were the most predominant, we wanted to determine the prevalence of the serogroups in cattle feces. Therefore, our secondary objectives were to, develop, validate, and use a mPCR assay to detect six predominant serogroups of non-top-7 STEC (O2, O74, O109, O131, O168 and O171) and utilize the assay to determine the prevalence of the predominant serogroups in feedlot cattle feces.

Materials and Methods

Non-top-7 STEC strains: A total of 359 STEC strains that were negative for the top-7 serogroups (O26, O45, O103, O11, O121, O145, and O157), obtained from two feedlot studies conducted in 2013 and 2014, were used. One-hundred-eighteen strains were isolated from twenty-four pen-floor fecal samples (n=576) collected weekly from two pens with 120-300 cattle each at a commercial feedlot in Nebraska during a 12-week period from June thru August of 2013 (Dewsbury et al., 2017). The remaining strains (n=241) were isolated from pen-floor fecal samples that were collected weekly from four-to-six pens with 30-450 cattle each at one commercial feedlot in Nebraska and one commercial feedlot in Texas during a 12-week period from June thru August of 2014 (n=1,886).

Both studies were designed to detect and estimate the prevalence of the top-7 STEC in cattle feces. For the detection and isolation of the top-7 STEC, serogroup-specific immunomagnetic beads were utilized to capture the top-7 STEC from the feces and then plated on selective media. Putative colonies from the selective media were subjected to a 10-plex PCR assay (Bai et al., 2012) that targeted seven serogroups and three major virulence genes (stx1, stx2, and eae). Isolates that were negative for the seven serogroups, but were positive one or both stx genes, were stored in Cryogenic beads (CryoCare™, Key Scientific Products, Round Rock, TX).

Detection of serogroups: Eleven sets of mPCR assays, previously developed and validated (Ludwig et al., 2017) to detect serogroups of STEC, other than the top-7, were utilized to detect the serogroups of the 359 strains. Strains frozen in cryogenic beads were streaked onto blood agar plates (Remel, Lenexa, KS) and incubated overnight at 37 C. Single colony of each strain

was picked, suspended in 50 µl distilled water, boiled for 10 min, centrifuged at 9,300 x g for 5 min, and the supernatant was used in mPCR assays.

Development and Validation of mPCR targeting O2, O74, O109, O131, O168, and O171

Primer design

Serogroup-specific wzx gene, which encodes for the transmembrane lipid transporter enzyme, or flippase, required for O-polysaccharide export, was targeted in this assay. The nucleotide sequence of the targeted gene for each of the six serogroups were obtained from NCBI GenBank, and a search was performed for additional sequences using BLASTN on NCBI's Nucleotide collection (nr/nt) database. Distinct amplicon sizes were used as targets for amplification by the primers. Primers were designed to match available O-serogroup sequences. The sequences for each serogroup were aligned using ClustalX version 2.0 (Larkin et al., 2007) and primers were designed to amplify the targets with distinct amplicon sizes that can be differentiated by capillary gel electrophoresis (Table 8).

PCR running conditions

Primer stocks were prepared in 1 x TE buffer (10mM Tris, brought to pH 8.0 with HCL with 1mM EDTA; Integrated DNA Technologies, Inc., Coralville, IA) at concentrations of 0.1 nM/ µl (stock solution). Concentrations and volumes of six-primer pairs were mixed, equally. The PCR running conditions consisted of 94°C initial denaturation for 5 min followed by 30 cycles of 94°C denaturation for 30 s, 72°C annealing for 30 s, 68°C extension for 75 s, and a final step of 68°C for 7 min. The final reaction mixture consisted of 20 µl and included 10 µl of IQ Multiplex Powermix (2x), 1µl of six pairs of primer mix (after mix, 8.3 pM/ul for each primer), 2 µl of

template, and 7 µl of water. The primer sequences and amplicon sizes for serogroups, O2, O74, O109, O131, O168, and O171 are shown in Table 3.

Validation of the mPCR assay with pure cultures and feces spiked with pure cultures

The applicability of the assay to detect serogroups in cattle feces was validated with fecal samples spiked with pure cultures of the six serogroups. Initially, five feedlot steer's fecal samples, obtained from pen-floor fecal pats, from the university Beef Cattle Research Center, were enriched (1 g in 9 ml, incubated at 40°C for 6 h) in *Escherichia coli* (EC) broth (Difco™, Becton, Dickinson Co., Sparks, MD). One milliliter of enriched fecal sample was boiled for 10 min and centrifuged at 9,300 x g for 5 mins. The DNA in the supernatant was purified with a GeneClean® Turbo Kit (MP Biomedicals, LLC, Solon, Ohio) and subjected to the 6-plex PCR assay. A fecal sample confirmed as negative for the six serogroups was used to spike with pure cultures in the validation assay. Six STEC strains (14.1652, 4558-1, 12662-2, 12.3205, 15.013, 1044-1 of O2, O74, O109, O131, O168, and O171, respectively), previously isolated from cattle feces, were used to spike fecal samples. The strains were grown in Luria-Bertani (LB) broth, and cell concentration (CFU/ml) was determined by spread-plate method. Serial ten-fold dilutions of each strain (10⁻⁰ to 10⁻⁸) were prepared in LB broth and 100 µl of 10⁻⁵, 10⁻⁶ and 10⁻⁷ dilutions were spread-plated onto blood agar plates (four plates per dilution) to determine initial cell concentrations (CFU/ml). Approximately 50 g of PCR-negative feces was suspended in 450 ml of EC broth and dispensed into sterile tubes (9.4 ml/tube). Serial dilutions (10⁻⁰ to 10⁻⁸) of pure cultures of each of the six serogroups were inoculated (100 µl) into 9.4 ml of EC broth. Spiked fecal samples were enriched (40°C for 6 h) in EC broth. Deoxyribonucleic acid was extracted from the enriched sample with a GeneClean® Turbo Kit, and subjected to the 6-plex PCR assay.

The experiment was repeated with a different fecal sample. The detection limit of the assay for *E. coli* serogroups O2, O74, O109, O131, O168, and O171 with DNA extracted directly from spiked cattle feces was $\sim 10^6$ CFU/g and $\sim 10^2$ CFU/g before and after enrichment in EC broth, respectively (Table 9).

Subtyping of *stx* genes

The subtyping of the *stx* genes of the 359 strains of non-top-7 STEC was according to Shridhar et al., 2017. A touchdown PCR method was used to amplify *stx* genes. The *stx1* amplification protocol included an initial denaturation at 94°C for 5 min, 10 cycles of touch-down PCR (denature: 94° C for 30 s, annealing: 56–51° C (Δ -0.5° C) for 30 s; and extension: 72°C for 1 min 45 s) followed by 30 cycles of regular PCR (denature: 94°C for 30 s, annealing: 51°C for 30 s; and extension: 72°C for 1 min 45 s). The *stx2* amplification protocol included an initial denaturation at 94°C for 5 min, 10 cycles of touch-down PCR (denature: 94° C for 30 s, annealing: 47–44° C (Δ -0.3°C) for 30 s; and extension: 72° C for 1 min 45 s) followed by 30 cycles of regular PCR (denature: 94°C for 30 s, annealing: 44°C for 30 s; and extension: 72°C for 1 min 45 s). Sequencing analysis was preformed using methods from Shridhar et al., 2017, and Scheutz et al., 2012. Products of PCR were visualized using a Qiaxcel capillary electrophoresis system (Qiagen, Valencia, CA) and purified using a QIAquick PCR purification kit (Qiagen). DNA concentration and purity were analyzed using a spectrophotometer (NanoDrop-Thermo Scientific, Wilmington, DE). Nucleotide sequencing was performed on PCR products and primers by Genewiz, Inc., (South Plainfield, NJ). Chromatogram data was visualized using the CLC Main Workbench software. Individual analysis of conflicts and secondary peaks were conducted for all sequences. Consensus sequences were produced. Intergenic sequences were removed, resulting in the translation of nucleotide sequences to amino acid sequences. Shiga

toxin subtypes were determined based on the amino acid motifs that define each *stx* subtype (Scheutz et al., 2012).

Results

In the first feedlot study, of the 118 strains tested, 95 (80.5%) belonged to 11 serogroups and 23 (19.5%) were unidentifiable (Table 1). The 11 identified serogroups included the following: O168 (29.7%), O109 (20.3%), O171 (9.3%), O74 (6.8%), O104 (4.2%), O2 (4.2%), O8 (2.5 %), O175 (.8%), O131 (0.8%), O160 (0.8%), and O169 (0.8%). In the second feedlot study, 203 strains (84.2%) of the 241 belonged to 19 serogroups and 38 (15.8%) were unidentified. The 19 serogroups included O168, O109, O131, O2, O104, O8, O74, O171, O136, O178, O102, O98, O175, O118, O113, O96, O76, O35, and O20. (Table 1). The seven strains positive for the *stx1* and *eae* genes were O118 (n=1), O98 (n=1), and O74 (n=5). The 36 strains positive for the *stx2* and *eae* genes were O175 (n=1) and O109 (n=35). A total of 16 strains belonging to O76 (n=1), O104 (n=12), and O109 (n=3) were positive for *stx1*. Out of the 241 strains, 162 strains belonging to O20 (n=1), O35 (n=1), O96 (n=1), O102 (n=2), O136 (n=3), O171 (n=4), O8 (n=1), O2 (n=11), O131 (n=28), O169 (n=72) and unknown (n=38) were positive for *stx2*. The 15 strains positive for both *stx1* and *stx2* were O113 (n=1), O178 (n=2), O8 (n=4), O2 (n=7), and O109 (n=1). (Table 7)

Detection of STEC strains from the two feedlot studies indicated that serogroups O2, O74, O109, O131, O168, and O171 were the predominant serogroups (Tables 6 and 7). In order to detect the prevalence of the six serogroups of non-top-7, a six-plex PCR assay was designed. The assay was validated with pure cultures and with feces spiked with pure cultures. The viable cell concentrations of strains of serogroups O2, O74, O109, O131, O168, and O171, cultured individually in LB broth, were 3.2×10^8 , 5.6×10^8 , 3.1×10^8 , 4.1×10^8 , 4.5×10^8 , and 6.5×10^8 ,

respectively. The detection limit of the assay for *E. coli* serogroups O2, O74, O109, O131, O168, and O171 with DNA extracted directly from spiked cattle feces was $\sim 10^6$ CFU/g and $\sim 10^2$ CFU/g before and after enrichment in EC broth, respectively (Table 9). A total of 911 fecal samples collected from commercial feedlots were subjected the six-plex PCR assay to determine the prevalence of the six predominant serogroups of non-top-7 STEC. Among the six serogroups of the major non-top 7 STEC, serogroups of O2 (59.4; 38.0), O109 (91.3; 84.9), O168 (78.4; 57.3), and O171 (86.8; 73.9) were the dominant serogroups in cattle feces collected from commercial feedlots from 2013 and 2014 studies, respectively (Tables 10 and 11). Only 29 (3.2%) fecal samples were negative for any of the six serogroups (Table 7). A majority of the fecal samples from the 2013 study (67.8%) were positive for 3 or 4 serogroups, and a majority from the 2014 study (51.5%) were positive for 2 or 3 serogroups (Table 12). Fecal samples that contained the six serogroups were also positive for the three major virulence genes, *stx1*, *stx2*, and *eae* (Tables 10 and 11).

Discussion

There are several pathotypes of intestinal pathogenic *E. coli*: Shiga toxin-producing *E. coli* (STEC)/ enterohemorrhagic *E. coli* (EHEC), enterotoxigenic *E. coli* (ETEC), enteropathogenic *E. coli* (EPEC), enteroinvasive *E. coli* (EIEC), enteroaggregative *E. coli* (EAEC), and diffusely adherent *E. coli* (DAEC). Both STEC and EHEC cause mild to severe hemorrhagic colitis with complications of HUS, particularly in children. Serotype O157:H7 is the most common STEC involved in foodborne infections, although other serogroups can also cause hemorrhagic colitis, including HUS. According to the National Enteric Disease Surveillance of STEC's 2014 annual report (updated as of Feb. 8th, 2017), incidence of

infection due to non-O157 serogroups continued to rise to 0.79 cases per 100,000 people whereas O157 incidence of infection continued to decrease to 0.60 cases per 100,000 people, in 2014.. Additionally, infections in infants (< 1 year of age) were mostly due to non-O157 serogroups (1.19 vs. 1.97 cases per 100,000 population for girls and 0.69 vs. 2.18 cases per 100,000 population for boys).

Shiga toxins are major virulence factors of STEC. While presence of serotype serves as an initial risk indicator, determining the Shiga toxin subtype of STEC strains provides a more reliable estimate of risk of human illness. Studies show that *stx2a* and *stx2c* subtypes represent most HUS cases in humans (Friedrich et al., 2002; Persson et al., 2007). Shiga toxin 2 subtypes, d, e, and f, have each been recognized as a severity indicator for illness in humans (Bielaszewska et al., 2006).

Escherichia coli hemolysin, also called α -hemolysin (Cavalieri et al., 1984), encoded by *ehxA*, is a major virulence factor that belongs to the repeats-in-toxin (RTX) family and lyses red blood cells in the host. The *hlyA* gene, also known as *ehxA*, encodes the toxin and is associated with *hlyC* and *hlyB/D* genes which encode for posttranslational modification, and components of the extracellular secretion system, respectively (Welch et al., 2005). Although much is unknown about the pathogenesis of Ehx/Hly, it is thought that iron from hemoglobin, released by mechanisms attributed to Hly, could stimulate the growth of STEC in the intestinal tract (Law et al., 1995).

In this study, prevalence of the major non top-7 STEC (O2, O74, O109, O131, O168, and O171) in commercial feedlot cattle was determined by developing and validating an end-point mPCR assay. Additionally, Shiga toxin subtypes associated with the major non top-7 STEC and others were identified based on restriction patterns after digestion of nucleotide sequences.

Among the six serogroups of non-top 7 STEC, serogroups of O2 (59.4; 38.0), O109 (91.3; 84.9), O168 (78.4; 57.3), and O171 (86.8; 73.9) were dominant in 2013 (%), and 2014 (%) studies, respectively. The most common subtypes carried by the dominant non top-7 STEC in this study were *stx2a* and *stx2c* for the 2013 and 2014 study which, as mentioned above, represent most HUS cases in humans; We found dominant serogroups O2, O109, O168, and O171 to be positive for *stx2a* at 72.0%, 56.3%, 87.9%, and 73.3%, respectively. Friedrich et al., (2002) found that one third of non-O157 serotypes positive for *stx2c* were also positive for *eae*, that 35.7% of HUS patients were positive for *stx2c*, and that the risk of developing HUS, post-infection, from the *stx2c* genotype was lower than the *stx2* genotype in STEC.

Isolating STEC from fecal samples in an efficient and effective manner is an important step in the overall process of evaluating risk of illness in humans, especially regarding HUS. Screening fecal samples by PCR for *stx*, followed by rapid subtyping of *stx2* genes in the presence or absence of *eae* provides physicians with trustworthy predictive values to assess the risk of their patients developing HUS from illness induced by STEC. This also provides a method of evaluation for those concerned about STEC shedding in cattle feces, which is a source of food and water contamination.

Table 6. Serogroup and virulence gene identification of enterohemorrhagic *Escherichia coli* and Shiga toxin-producing *Escherichia coli* strains isolated from commercial feedlot cattle feces (2013 study) by end-point PCR.

Serogroup	Total positive (%)	Enterohemorrhagic <i>E. coli</i>				Shiga toxin-producing <i>E. coli</i>		
		<i>stx1</i>	<i>stx2</i>	<i>eae</i>	<i>ehxA</i>	<i>stx1</i>	<i>stx2</i>	<i>ehxA</i>
O131	1 (0.8)	-	-	-	-	-	1 <i>stx2d</i>	1
O160	1 (0.8)	-	-	-	-	-	1 <i>stx2a</i>	-
O169	1 (0.8)	-	-	-	-	-	1 <i>stx2a</i>	1
O175	1 (0.8)	-	-	-	-	-	1 <i>stx2d</i>	1
O8	3 (2.5)	-	-	-	1	1 <i>stx1a</i>	2 <i>stx2a</i> ; 1 <i>stx2d</i>	1
O2	5 (4.2)	-	-	-	-	1 <i>stx1a</i>	4 <i>stx2a</i> ; 1 <i>stx2c</i>	5
O104	5 (4.2)	-	-	-	-	5 <i>stx1c</i>	-	5
O74	8 (6.8)	8 <i>stx1a</i>	-	8	8	-	-	-
O171	11 (9.3)	1 <i>stx1a</i>	-	1	-	-	9 <i>stx2a</i> ; 1 <i>stx2c</i>	2
	24		12 <i>stx2a</i> ;					
O109	(20.3)	-	4 <i>stx2c</i>	16	16	8 <i>stx1a</i>	-	-
	35						31 <i>stx2a</i> ; 1	
O168	(29.7)	-	-	-	-	1 <i>stx1a</i>	<i>stx2c</i> ; 2 <i>stx2d</i>	20
	23						16 <i>stx2a</i> ; 2	
Unknown	(19.5)	-	-	-	1	1 <i>stx1a</i>	<i>stx2c</i> ; 4 <i>stx2d</i>	22

Table 7. Serogroup and virulence gene identification of enterohemorrhagic *Escherichia coli* and Shiga toxin-producing *Escherichia coli* strains isolated from commercial feedlot cattle feces (2014 study) by end-point PCR.

Serogroup	Total positive (%)	Enterohemorrhagic <i>E. coli</i>				Shiga toxin-producing <i>E. coli</i>		
		<i>stx1</i>	<i>stx2</i>	<i>eae</i>	<i>ehxA</i>	<i>stx1</i>	<i>stx2</i>	<i>ehxA</i>
O20	1 (.4)	-	-	-	-	-	1 <i>stx2a</i>	1
O35	1 (.4)	-	-	-	-	-	1 <i>stx2a</i>	-
O76	1 (.4)	-	-	-	-	1 <i>stx1a</i>	-	-
O96	1 (.4)	-	-	-	-	-	1 <i>stx2a</i>	1
O113	1 (.4)	-	-	-	-	1 <i>stx1a</i>	1 <i>stx2a</i>	1
O118	1 (.4)	1 <i>stx1a</i>	-	1	1	-	-	-
O175	1 (.4)	-	1 <i>stx2a</i>	1	1	-	-	-
O98	1 (.4)	1 <i>stx1a</i>	-	1	1	-	-	-
O102	2 (.8)	-	-	-	-	-	2 <i>stx2a</i>	2
O178	2 (.8)	-	-	-	-	2 <i>stx1a</i>	2 <i>stx2a</i>	2
O136	3 (1.3)	-	-	-	-	-	2 <i>stx2d</i> ; 1 <i>stx2a</i>	3
O171	4 (1.7)	-	-	-	-	-	2 <i>stx2a</i> ; 1 <i>stx2c</i> ; 1 <i>stx2d</i>	-
O74	5 (2.1)	5 <i>stx1a</i>	-	5	5	-	-	-
O8	7 (2.9)	-	-	-	-	4 <i>stx1a</i>	4 <i>stx2a</i> ; 1 <i>stx2d</i>	5
O104	12 (5.0)	-	-	-	-	11 <i>stx1c</i> ; 1 <i>stx1a</i>	-	12
O2	20 (8.3)	-	-	-	-	7 <i>stx1a</i>	14 <i>stx2a</i> ; 1 <i>stx2c</i> ; 3 <i>stx2d</i>	16
O131	28 (11.6)	-	-	-	-	-	19 <i>stx2a</i> ; 9 <i>stx2c</i>	1

			24 <i>stx2a</i> ;					
	40	-	5 <i>stx2c</i> ; 7					-
O109	(16.6)		<i>stx2d</i>	36	35	3 <i>stx1a</i>	1 <i>stx2a</i>	
O168	72 (29.9)	-	-	-	-	-	63 <i>stx2a</i> ; 6 <i>stx2d</i> ; 3 <i>stx2c</i>	53
Unknown	37 (15.8)	-	-			-	27 <i>stx2a</i> ; 4 <i>stx2c</i> ; 6 <i>stx2d</i>	18

Table 8. Target genes, primers sequences, and amplicon sizes to detect six major non top-7 Shiga toxin-producing *Escherichia coli*.

Genes	Primer	Primers sequence	Amplicon size
wzxO2/O50	O2-F	TGGCCTTGTTTCGATATACTGCGGA	819
	O2-R	TCACGAGCTGAGCGAAACTGTTCA	
wzxO131	O131-F	GTGATTTCTGGGGCAACATT	655
	O131-R	AAGCCTGCCCTAAACAAAGC	
wzxO168	O168-F	TGTCGACTTTGGGAAATGTGG	336
	O168-R	CTGCAGAGGCCAATTCAGGT	
wzxO74	O74-F	CTGGTCAATGGCAAGCTGTA	303
	O74-R	ATGCAAAAATCCAAGCCAAT	
wzxO171	O171-F	TGCTCAAGTGGCATGCAGAT	281
	O171-R	TGCAACCTGATATCCAGCAGT	
wzxO109	O109-F	TCTCTCTCGACATACCCGCGCTT	204
	O109-R	ACCGTAGCCCAAAGAGCCACA	

Table 9. Fecal spike validation of six major non top-7 Shiga toxin-producing *Escherichia coli* from cattle feces by end-point PCR method of detection.

Concentrations		Serogroups				
Pre-Enrichment	O2	O74	O109	O131	O168	O171
10 ⁶	+	+	+	+	+	+
10 ⁵	+	-	+	+	+	+
10 ⁴	-	-	+	+	+	+
10 ³	-	-	+	-	+	+
10 ²	-	-	+	-	-	-
10 ¹	-	-	+	-	-	-
10 ⁰	-	-	-	-	-	-
Post-Enrichment						
10 ⁶	+	+	+	+	+	+
10 ⁵	+	+	+	+	+	+
10 ⁴	+	+	+	+	+	+
10 ³	+	-	+	+	+	+
10 ²	+	+	+	+	+	+
10 ¹	-	-	+	-	+	+
10 ⁰	-	-	+	-	+	+

Table 10. Detection of six major non top-7 Shiga toxin-producing *Escherichia coli* and virulence genes, *stx1*, *stx2*, and *eae*, from cattle feces (2013 study) in a commercial feedlot setting by end-point PCR.

Serogroup	Total positive (%)	<i>stx1</i>	<i>stx2</i>	<i>eae</i>
O2	341 (59.4)	245 (71.8)	326 (95.6)	331 (97.6)
O74	113 (19.7)	83 (73.5)	107 (94.7)	111 (98.2)
O109	524 (91.3)	346 (66.0)	501 (95.6)	512 (97.7)
O131	10 (1.7)	5 (50.0)	10 (100.0)	10 (100.0)
O168	450(78.4)	310 (68.9)	437 (97.1)	442 (98.2)
O171	498 (86.8)	332 (66.7)	478 (96.0)	488 (98.0)

Table 11. Prevalence of six major non top-7 Shiga toxin-producing *Escherichia coli* and virulence genes, *stx1*, *stx2*, and *eae*, from cattle feces (2014 study) in a commercial feedlot setting by end-point PCR.

Serogroup	Total Positive (%)	<i>stx1</i>	<i>stx2</i>	<i>eae</i>
O2	128 (38.0)	85 (66.4)	108 (84.4)	126 (98.4)
O74	56 (16.6)	38 (67.9)	41 (73.2)	51 (91.1)
O109	286 (84.9)	157 (54.9)	220 (76.9)	263 (92.0)
O131	8 (2.4)	7 (87.5)	8 (100.0)	8 (100.0)
O168	193 (57.3)	120 (62.2)	167 (86.5)	181 (93.8)
O171	249 (73.9)	146 (58.6)	199 (79.9)	240 (96.4)

Table 12. Number and percentage of samples testing positive to multiple major non top-7 Shiga toxin-producing *Escherichia coli* serogroups ((O2, O74, 0109, 0131, 0168, 0171) by end-point PCR.

Study (no. samples tested)	Total positive (%)	Number of serogroups (%)						
		0	1	2	3	4	5	6
2013 (n=574)	560	14	33	63	156	233	73	2
	(97.6)	(2.4)	(5.7)	(11)	(27.2)	(40.6)	(12.7)	(0.3)
2014 (n=337)	322	15	46	87	87	73	27	2
	(95.5)	(4.5)	(13.6)	(25.8)	(25.8)	(21.7)	(8.0)	(0.6)

References

- Aidar-Ugrinovich, L., Blanco, J., Blanco, M., Blanco, J. E., Leomil, L., Dahbi, G., & de Castro, A. P. (2007). Serotypes, virulence genes, and intimin types of Shiga toxin-producing *Escherichia coli* (STEC) and enteropathogenic *E. coli* (EPEC) isolated from calves in São Paulo, Brazil. *International Journal of Food Microbiology*, 115(3), 297-306.
- Bai, Justin Ludwig; X. Shi; K. Capps; E. Roberts; C. DebRoy; R. Phebus; T.G. Nagaraja; and J. Development and validation of multiplex PCR assays to identify 'Non-Top-7' serogroups of Shiga toxin-producing *Escherichia coli* known to be present in cattle feces. (2017).
- Bai, J., Paddock, Z. D., Shi, X., Li, S., An, B., & Nagaraja, T. G. (2012). Applicability of a multiplex PCR to detect the seven major Shiga toxin-producing *Escherichia coli* based on genes that code for serogroup-specific O-antigens and major virulence factors in cattle feces. *Foodborne Pathogens and Disease*, 9(6), 541-548.
- Basu, D., X. P. Li, J. N. Kahn, K. L. May, P. C. Kahn, and N. E. Tumer. The A1 subunit of Shiga toxin 2 has higher affinity for ribosomes and higher catalytic activity than the A1 Subunit of Shiga toxin 1. *Infect Immun* 84, no. 1 (Oct 19 2015): 149-61.
- Bettelheim, K. A. (2007). The non-O157 Shiga-toxigenic (verocytotoxigenic) *Escherichia coli*; under-rated pathogens. *Critical Reviews In Microbiology*, 33(1), 67-87.
- Beutin, L., Miko, A., Krause, G., Pries, K., Haby, S., Steege, K., & Albrecht, N. (2007). Identification of human-pathogenic strains of Shiga toxin-producing *Escherichia coli* from food by a combination of serotyping and molecular typing of Shiga toxin genes. *Applied And Environmental Microbiology*, 73(15), 4769-4775.
- Bhakdi, S., M. Muhly, S. Korom, and G. Schmidt. Effects of *Escherichia coli* hemolysin on human monocytes. cytotoxic action and stimulation of interleukin 1 release. *J Clin Invest* 85,: 1746-53.
- Bielaszewska, M., A. W. Friedrich, T. Aldick, R. Schurk-Bulgrin, and H. Karch. Shiga toxin activatable by intestinal mucus in *Escherichia coli* isolated from humans: predictor for a severe clinical outcome. *Clin Infect Dis* 43,: 1160-7.
- Blanco, M., Blanco, J. E., Mora, A., Dahbi, G., Alonso, M. P., González, E. A., ... & Blanco, J. (2004). Serotypes, virulence genes, and intimin types of Shiga toxin (verotoxin)-producing *Escherichia coli* isolates from cattle in Spain and identification of a new intimin variant gene (eae-ξ). *Journal of Clinical Microbiology*, 42(2), 645-651.
- Blanco, M., Padola, N. L., Krüger, A., Sanz, M. E., Blanco, J. E., González, E. A., & Arroyo, G. H. (2004). Virulence genes and intimin types of Shiga-toxin-producing *Escherichia coli* isolated from cattle and beef products in Argentina. *International Microbiology*, 7(4), 269-276.

- Brandal, L. T., A. L. Wester, H. Lange, I. Lobersli, B. A. Lindstedt, L. Vold, and G. Kapperud. Shiga toxin-producing *Escherichia coli* infections in Norway, 1992-2012: characterization of isolates and identification of risk factors for haemolytic uremic syndrome. *BMC Infect Dis* 15: 324.
- Brooks, J. T., E. G. Sowers, J. G. Wells, K. D. Greene, P. M. Griffin, R. M. Hoekstra, and N. A. Strockbine. Non-O157 Shiga toxin-producing *Escherichia coli* infections in the United States, 1983-2002. *J Infect Dis* 192, : 1422-9.
- Cavaleri, S. J., G. A. Bohach, and I. S. Snyder. *Escherichia coli* alpha-hemolysin: characteristics and probable role in pathogenicity. *Microbiol Rev* 48,: 326-43.
- Cergole-Novella, M. C., Nishimura, L. S., Dos Santos, L. F., Irino, K., Vaz, T. M. I., Bergamini, A. M., & Guth, B. E. C. (2007). Distribution of virulence profiles related to new toxins and putative adhesins in Shiga toxin-producing *Escherichia coli* isolated from diverse sources in Brazil. *FEMS Microbiology Letters*, 274(2), 329-334.
- Coimbra, R. S., Grimont, F., Lenormand, P., Burguière, P., Beutin, L., & Grimont, P. A. (2000). Identification of *Escherichia coli* O-serogroups by restriction of the amplified O-antigen gene cluster (rfb-RFLP). *Research in Microbiology*, 151(8), 639-654.
- Cull, C. A., D. G. Renter, D. M. Dewsbury, L. W. Noll, P. B. Shridhar, S. E. Ives, T. G. Nagaraja, and N. Cernicchiaro. Feedlot- and pen-level prevalence of Enterohemorrhagic *Escherichia coli* in feces of commercial feedlot cattle in two major U.S. cattle feeding areas. *Foodborne Pathogens Disease*. 14, : 309-17.
- DebRoy, C., Fratamico, P. M., Yan, X., Baranzoni, G., Liu, Y., Needleman, D. S., & Mwangi, M. (2016). Comparison of O-antigen gene clusters of all O-serogroups of *Escherichia coli* and proposal for adopting a new nomenclature for O-typing. *PLoS One*, 11(1), e0147434.
- Dewsbury, D. M., Renter, D. G., Shridhar, P. B., Noll, L. W., Shi, X., Nagaraja, T. G., & Cernicchiaro, N. (2015). Summer and winter prevalence of Shiga toxin-producing *Escherichia coli* (STEC) O26, O45, O103, O111, O121, O145, and O157 in feces of feedlot cattle. *Foodborne Pathogens and Disease*, 12(8), 726-732.
- Ennis, C., McDowell, D., & Bolton, D. J. (2012). The prevalence, distribution and characterization of Shiga toxin-producing *Escherichia coli* (STEC) serotypes and virulotypes from a cluster of bovine farms. *Journal of Applied Microbiology*, 113(5), 1238-1248.
- Ethelberg, S., K. E. Olsen, F. Scheutz, C. Jensen, P. Schiellerup, J. Enberg, A. M. Petersen, *et al.* "virulence factors for hemolytic uremic syndrome, Denmark." *Emerg Infect Dis* 10, : 842-7.
- Fagan, P. K., Hornitzky, M. A., Bettelheim, K. A., & Djordjevic, S. P. (1999). Detection of Shiga-like toxin (*stx1* and *stx2*), intimin (*eaeA*), and enterohemorrhagic *Escherichia coli* (EHEC)

- hemolysin (EHEC *hlyA*) genes in animal feces by multiplex PCR. *Applied and Environmental Microbiology*, 65(2), 868-872.
- Feng, P. C., and S. Reddy. Prevalences of Shiga toxin subtypes and selected other virulence factors among Shiga-toxigenic *Escherichia coli* strains isolated from fresh produce." *Appl Environ Microbiol* 79, : 6917-23.
- Fratamico, P. M., & Bagi, L. K. (2012). Detection of Shiga toxin-producing *Escherichia coli* in ground beef using the GeneDisc real-time PCR system. *Frontiers in Cellular and Infection Microbiology*, 2.
- Fratamico, P. M., Bagi, L. K., Cray Jr, W. C., Narang, N., Yan, X., Medina, M., & Liu, Y. (2011). Detection by multiplex real-time polymerase chain reaction assays and isolation of Shiga toxin-producing *Escherichia coli* serogroups O26, O45, O103, O111, O121, and O145 in ground beef. *Foodborne Pathogens and Disease*, 8(5), 601-607.
- Fratamico, P. M., Wasilenko, J. L., Garman, B., DeMARCO, D. R., Varkey, S., Jensen, M., & Tice, G. (2014). Evaluation of a multiplex real-time PCR method for detecting Shiga toxin-producing *Escherichia coli* in beef and comparison to the US Department of Agriculture Food Safety and Inspection Service microbiology laboratory guidebook method. *Journal of Food Protection*, 77(2), 180-188.
- Friedrich, A. W., M. Bielaszewska, W. L. Zhang, M. Pulz, T. Kuczius, A. Ammon, and H. Karch. *Escherichia coli* harboring Shiga toxin 2 gene variants: frequency and association with clinical symptoms. *J Infect Dis* 185, : 74-84.
- Fukushima, H., & Seki, R. (2004). High numbers of Shiga toxin-producing *Escherichia coli* found in bovine faeces collected at slaughter in Japan. *FEMS Microbiology Letters*, 238(1), 189-197.
- Gannon, V. P., Rashed, M., King, R. K., & Thomas, E. J. (1993). Detection and characterization of the *eae* gene of Shiga-like toxin-producing *Escherichia coli* using polymerase chain reaction. *Journal of Clinical Microbiology*, 31(5), 1268-1274.
- Gonzalez, A. G., Cerqueira, A. M., Guth, B. E., Coutinho, C. A., Liberal, M. H. T., Souza, R. M., & Andrade, J. R. (2016). Serotypes, virulence markers and cell invasion ability of Shiga toxin-producing *Escherichia coli* strains isolated from healthy dairy cattle. *Journal Of Applied Microbiology*, 121(4), 1130-1143.
- Guth, B. E., Chinen, I., Miliwebsky, E., Cerqueira, A. M., Chillemi, G., Andrade, J. R., & Rivas, M. (2003). Serotypes and Shiga toxin genotypes among *Escherichia coli* isolated from animals and food in Argentina and Brazil. *Veterinary Microbiology*, 92(4), 335-349.
- Hara-Kudo, Y., Konishi, N., Ohtsuka, K., Iwabuchi, K., Kikuchi, R., Isobe, J., ... & Tanouchi, A. (2016). An interlaboratory study on efficient detection of Shiga toxin-producing *Escherichia coli* O26,

- O103, O111, O121, O145, and O157 in food using real-time PCR assay and chromogenic agar. *International Journal of Food Microbiology*, 230, 81-88.
- Hornitzky, M. A., Bettelheim, K. A., & Djordjevic, S. P. (2001). The detection of Shiga toxin-producing *Escherichia coli* in diagnostic bovine faecal samples using vancomycin-cefixime-cefsulodin blood agar and PCR. *FEMS Microbiology Letters*, 198(1), 17-22.
- Hornitzky, M. A., Mercieca, K., Bettelheim, K. A., & Djordjevic, S. P. (2005). Bovine feces from animals with gastrointestinal infections are a source of serologically diverse atypical enteropathogenic *Escherichia coli* and Shiga toxin-producing *E. coli* strains that commonly possess intimin. *Applied and Environmental Microbiology*, 71(7), 3405-3412.
- Hughes, J. M., Wilson, M. E., Johnson, K. E., Thorpe, C. M., & Sears, C. L. (2006). The emerging clinical importance of non-O157 Shiga toxin—producing *Escherichia coli*. *Clinical Infectious Diseases*, 43(12), 1587-1595.
- Hussein, H. S., & Bollinger, L. M. (2005). Prevalence of Shiga toxin-producing *Escherichia coli* in beef. *Meat Science*, 71(4), 676-689.
- Iguchi, A., Iyoda, S., Kikuchi, T., Ogura, Y., Katsura, K., Ohnishi, M., ... & Thomson, N. R. (2014). A complete view of the genetic diversity of the *Escherichia coli* O-antigen biosynthesis gene cluster. *DNA Research*, 22(1), 101-107.
- Iguchi, A., Iyoda, S., Seto, K., Morita-Ishihara, T., Scheutz, F., & Ohnishi, M. (2015). *Escherichia coli* O-genotyping PCR: a comprehensive and practical platform for molecular O serogrouping. *Journal of Clinical Microbiology*, 53(8), 2427-2432.
- Irino, K., Kato, M. A. M. F., Vaz, T. M. I., Ramos, I. I., Souza, M. A. C., Cruz, A. S., & Guth, B. E. C. (2005). Serotypes and virulence markers of Shiga toxin-producing *Escherichia coli* (STEC) isolated from dairy cattle in São Paulo State, Brazil. *Veterinary Microbiology*, 105(1), 29-36.
- Joensen, K. G., Tetzschner, A. M., Iguchi, A., Aarestrup, F. M., & Scheutz, F. (2015). Rapid and easy in silico serotyping of *Escherichia coli* using whole genome sequencing (WGS) data. *Journal of Clinical Microbiology*, .00008.
- Jurgens, D., M. Ozel, and N. B. Takaisi-Kikuni. Production and characterization of *Escherichia coli* enterohemolysin and its effects on the structure of erythrocyte membranes." *Cell Biol Int* 26, : 175-86.
- Karch, H., E. Denamur, U. Dobrindt, B. B. Finlay, R. Hengge, L. Johannes, E. Z. Ron, *et al.* "The enemy within us: lessons from the 2011 european *Escherichia coli* O104:H4 outbreak." *EMBO Mol Med* 4, : 841-8.
- Karmali, M. A. (2004). Prospects for preventing serious systemic toxemic complications of Shiga toxin—producing *Escherichia coli* infections using Shiga toxin receptor analogues. *The Journal of Infectious Diseases*, 189(3), 355-359.

- Keane, W. F., R. Welch, G. Gekker, and P. K. Peterson. Mechanism of *Escherichia coli* alpha-hemolysin-induced injury to isolated renal tubular cells. *Am J Pathol* 126,: 350-7.
- Kido, N., & Kobayashi, H. (2000). A single amino acid substitution in a mannosyltransferase, WbdA, converts the *Escherichia coli* O9 polysaccharide into O9a: generation of a new O-serotype group. *Journal of Bacteriology*, 182(9), 2567-2573.
- Knirel, Y. A., Qian, C., Shashkov, A. S., Sizova, O. V., Zdorovenko, E. L., Naumenko, O. I., & Liu, B. (2016). Structural relationships between genetically closely related O-antigens of *Escherichia coli* and *Shigella* spp. *Biochemistry*, 81(6), 600.
- Lacher, D. W., Gangiredla, J., Jackson, S. A., Elkins, C. A., & Feng, P. C. (2014). Novel microarray design for molecular serotyping of Shiga toxin-producing *Escherichia coli* strains isolated from fresh produce. *Applied and Environmental Microbiology*, 80(15), 4677-4682.
- Larkin, M. A., Blackshields, G., Brown, N. P., Chenna, R., McGettigan, P. A., McWilliam, H., & Thompson, J. D. (2007). Clustal W and Clustal X version 2.0. *Bioinformatics*, 23(21), 2947-2948.
- Law, D., and J. Kelly. Use of heme and hemoglobin by *Escherichia coli* O157 and other Shiga-like-toxin-producing *E. coli* serogroups. *Infect Immun* 63, : 700-2.
- Li, H., A. Granat, V. Stewart, and J. R. Gillespie. Rpos, H-Ns, and Dsra Influence Ehec Hemolysin Operon (Ehxcabd) Transcription in *Escherichia coli* O157:H7 Strain Edl933." *FEMS Microbiol Lett* 285, : 257-62.
- Liu, Y., & Fratamico, P. (2006). *Escherichia coli* O antigen typing using DNA microarrays. *Molecular and Cellular Probes*, 20(3), 239-244.
- Liu, Y., Yan, X., DebRoy, C., Fratamico, P. M., Needleman, D. S., Li, R. W., ... & Toro, M. (2015). *Escherichia coli* O-antigen gene clusters of serogroups O62, O68, O131, O140, O142, and O163: DNA sequences and similarity between O62 and O68, and PCR-based serogrouping. *Biosensors*, 5(1), 51-68.
- Meichtri, L., Miliwebsky, E., Gioffré, A., Chinen, I., Baschkier, A., Chillemi, G., & Rivas, M. (2004). Shiga toxin-producing *Escherichia coli* in healthy young beef steers from Argentina: prevalence and virulence properties. *International Journal Of Food Microbiology*, 96(2), 189-198.
- Miko, A., Rivas, M., Bentancor, A., Delannoy, S., Fach, P., & Beutin, L. (2014). Emerging types of Shiga toxin-producing *E. coli* (STEC) O178 present in cattle, deer, and humans from Argentina and Germany. *Frontiers in Cellular and Infection Microbiology*, 4.

- Monaghan, A., Byrne, B., Fanning, S., Sweeney, T., McDowell, D., & Bolton, D. J. (2011). Serotypes and virulence profiles of non-O157 Shiga toxin-producing *Escherichia coli* isolates from bovine farms. *Applied and Environmental Microbiology*, 77(24), 8662-8668.
- Noll, L. W., P. B. Shridhar, D. M. Dewsbury, X. Shi, N. Cernicchiaro, D. G. Renter, and T. G. Nagaraja (2015). "A comparison of culture- and PCR-based methods to detect six major non-O157 serogroups of Shiga toxin-producing *Escherichia coli* in cattle feces." *PLoS One* 10, : e0135446.
- Noll, L. W., P. B. Shridhar, X. Shi, B. An, N. Cernicchiaro, D. G. Renter, T. G. Nagaraja, and J. Bai (2015). A four-plex real-Time PCR assay, based on *rfbe*, *stx1*, *stx2*, and *eae* genes, for the detection and quantification of Shiga toxin-producing *Escherichia coli* O157 in cattle feces. *Foodborne Pathog Dis* 12, : 787-94.
- Norman, K. N., Clawson, M. L., Strockbine, N. A., Mandrell, R. E., Johnson, R., Ziebell, K., & Bono, J. L. (2015). Comparison of whole genome sequences from human and non-human *Escherichia coli* O26 strains. *Frontiers in Cellular and Infection Microbiology*, 5.
- Ogura, Y., T. Ooka, A. Iguchi, H. Toh, M. Asadulghani, K. Oshima, T. Kodama, *et al.* (2009). Comparative Genomics reveal the mechanism of the parallel evolution of O157 and non-O157 enterohemorrhagic *Escherichia coli*. *Proc Natl Acad Sci U S A* 106): 17939-44.
- Oliveira, M. G., Brito, J. R. F., Gomes, T. A. T., Guth, B. E. C., Vieira, M. A. M., Naves, Z. V. F. & Irino, K. (2008). Diversity of virulence profiles of Shiga toxin-producing *Escherichia coli* serotypes in food-producing animals in Brazil. *International Journal of Food Microbiology*, 127(1), 139-146.
- Padola, N. L., Sanz, M. E., Blanco, J. E., Blanco, M., Blanco, J., Etcheverria, A. I., ... & Parma, A. E. (2004). Serotypes and virulence genes of bovine Shiga toxigenic *Escherichia coli* (STEC) isolated from a feedlot in Argentina. *Veterinary Microbiology*, 100(1), 3-9.
- Persson, S., K. E. Olsen, S. Ethelberg, and F. Scheutz. Subtyping method for *Escherichia coli* Shiga toxin (Verocytotoxin) 2 variants and correlations to clinical manifestations. *J Clin Microbiol* 45, : 2020-4.
- Prager, R., A. Fruth, U. Siewert, U. Strutz, and H. Tschape. *Escherichia coli* encoding Shiga toxin 2f as an emerging human pathogen. *Int J Med Microbiol* 299, : 343-53.
- Saitoh, T., S. Iyoda, S. Yamamoto, Y. Lu, K. Shimuta, M. Ohnishi, J. Terajima, and H. Watanabe. Transcription of the *ehx* enterohemolysin gene is positively regulated by *grla*, a global regulator encoded within the locus of enterocyte effacement in enterohemorrhagic *Escherichia coli*. *J Bacteriol* 190, : 4822-30.
- Scheffer, J., W. Konig, J. Hacker, and W. Goebel. "Bacterial adherence and hemolysin production from *Escherichia coli* induces histamine and leukotriene release from various cells." *Infect Immun* 50,: 271-8.

- Scott, L., McGee, P., Walsh, C., Fanning, S., Sweeney, T., Blanco, J., & Sheridan, J. J. (2009). Detection of numerous verotoxigenic *E. coli* serotypes, with multiple antibiotic resistance from cattle faeces and soil. *Veterinary Microbiology*, 134(3), 288-293.
- Shridhar, P. B., C. Siepker, L. W. Noll, X. Shi, T. G. Nagaraja, and J. Bai. (2017). Shiga toxin subtypes of non-O157 *Escherichia coli* serogroups isolated from cattle feces. *Front Cell Infect Microbiol* 7: 121.
- Shridhar, P. B., L. W. Noll, X. Shi, N. Cernicchiaro, D. G. Renter, J. Bai, and T. G. Nagaraja. (2016). *Escherichia coli* O104 in feedlot cattle feces: prevalence, isolation and characterization. *PLoS One* 11, no. 3: e0152101.
- Slanec, T., Fruth, A., Creuzburg, K., & Schmidt, H. (2009). Molecular analysis of virulence profiles and Shiga toxin genes in food-borne Shiga toxin-producing *Escherichia coli*. *Applied and Environmental Microbiology*, 75(19), 6187-6197.
- Stenutz, R., Weintraub, A., & Widmalm, G. (2006). The structures of *Escherichia coli* O-polysaccharide antigens. *FEMS Microbiology Reviews*, 30(3), 382-403.
- Strockbine, N. A., L. R. Marques, J. W. Newland, H. W. Smith, R. K. Holmes, and A. D. O'Brien. Two toxin-converting phages from *Escherichia coli* O157:H7 strain 933 encode antigenically distinct toxins with similar biologic activities. *Infect Immun* 53, : 135-40.
- Suttorp, N., B. Floer, H. Schnittler, W. Seeger, and S. Bhakdi. Effects of *Escherichia coli* hemolysin on endothelial cell function. *Infect Immun* 58, : 3796-801.
- Tarr, P. I., Gordon, C. A., & Chandler, W. L. (2005). Shiga-toxin-producing *Escherichia coli* and haemolytic uraemic syndrome. *The Lancet*, 365(9464), 1073-1086.
- Tozzi, A. E., Caprioli, A., Minelli, F., Gianviti, A., De Petris, L., Edefonti, A., ... & Rizzoni, G. (2003). Shiga toxin-producing *Escherichia coli* infections associated with hemolytic uremic syndrome, Italy, 1988–2000. *Emerging Infectious Diseases*, 9(1), 106.
- Wasilenko, J. L., Frataamico, P. M., Sommers, C., DeMarco, D. R., Varkey, S., Rhoden, K., & Tice, G. (2014). Detection of Shiga toxin-producing *Escherichia coli* (STEC) O157: H7, O26, O45, O103, O111, O121, and O145, and *Salmonella* in retail raw ground beef using the DuPont™ BAX® system. *Frontiers in Cellular and Infection Microbiology*, 4.
- Weinstein, D. L., M. P. Jackson, J. E. Samuel, R. K. Holmes, and A. D. O'Brien. Cloning and sequencing of a Shiga-like toxin type II variant from *Escherichia coli* strain responsible for rdema disease of swine. *J Bacteriol* 170, : 4223-30.
- Welch, R. A. The *Escherichia coli* Hemolysin. *EcoSal Plus* 1, no. 2 (Nov 2005). "Pore-forming cytolysins of gram-negative bacteria. *Mol Microbiol* 5, no. 3 (Mar 1991): 521-8.

- Werber, D., Beutin, L., Pichner, R., Stark, K., & Fruth, A. (2008). Shiga toxin–producing *Escherichia coli* serogroups in food and patients, Germany. *Emerging Infectious Diseases*, 14(11), 1803.
- Woods, T. A., H. M. Mendez, S. Ortega, X. Shi, D. Marx, J. Bai, R. A. Moxley, *et al.* (2016). Development of 11-plex mol-pcr assay for the rapid screening of samples for Shiga toxin-producing *Escherichia coli*. *Front Cell Infect Microbiol* 6: 92.
- Zweifel, C., Schumacher, S., Blanco, M., Blanco, J. E., Tasara, T., Blanco, J., & Stephan, R. (2005). Phenotypic and genotypic characteristics of non-O157 Shiga toxin-producing *Escherichia coli* (STEC) from Swiss cattle. *Veterinary Microbiology*, 105(1), 37-45.